


INVACARE®
COMPASS™ SPT
TILT-IN-SPACE WHEELCHAIR



The Invacare® Compass™ Series adult manual wheelchairs provide the optimal solution for complex positioning and seating situations. The **Invacare® Compass™ SPT Wheelchair** combines a traditional tilt-in-space design with self-propelling functionality, all in one wheelchair.

Features

Pin and saddle style front riggings offered to be compatible with other self-propelling manual wheelchairs

Exclusive Invacare® Quick Set™ washers snap onto frame anywhere frequent adjustments might occur

Discrete tilt trigger lock-out levers help prevent people from accidentally disengaging the tilt mechanisms


Fold-down back and quick-release wheel option allow for easier transportability

Designed adjustable angle footplates have a dual attachment point for increased durability and a positive locking mechanism for dorsi and planar flexion adjustment

Adjustable angle stroller handle option offered for easier dependent propulsion

Depth adjustable frame – three different seat rail lengths available with at least 3" of adjustment on each frame

HCPCS Code: E1161



1. Front wheel configuration
2. Adjustable angle footplate
3. Dual gas cylinders
4. Quick set washer
5. Tilt trigger lock-out lever



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 12-250
120365

Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Tilt Range	Back Angle	Front Frame Angle	Tubing Diameter
SPT	Seat width + 9.5"	14"-22"W 15"-22"D	14.5"- 18.5" 1" increments	0°-45°	90°-120° 10° increments	N/A	1"
	Knee-to-Heel Range	Wheels	Casters	Total Weight*	Weight Capacity	Limited Warranty	
	13.5"-21.5" 1" increments	20", 22", 24"	6", 8"	35 lb.	250 lb.	5 Years Frame/ Crossbraces	

* 16"x16" Base chair without accessories.

Please see the Owner's Manual or Invacare Configurator for complete option/specification combinations.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



Yes, you can.®