



twion[®]
ENERGIZE YOUR WHEELCHAIR

alber[®]

Simply mobile.

MADE IN GERMANY 





PUSHING THE LIMITS.

In today's world, one thing matters more than anything: mobility. That's why we need a mode of transport that is adapted to our dynamic lifestyles: efficient, smart and fun to drive.

Innovative drive technology from Alber delivers the efficient and flexible solution for future mobility.

The new twion power assist offers new freedom and maximum independence. The electric motors equipped with leading-edge digital electronics provide you with extra power for every propelling movement. The sensor fitted to the hand rim calculates precisely the right amount of assist.

The silent direct drive ensures high agility and simple handling. You are actively on the move and reach your destination quickly, efficiently and calmly.



FEEL THE SPEED.

Faster and further: The dynamic motors accelerate the wheelchair up to 6 mph. So you can really get going and enjoy the ride. The efficient energy recovery system feeds energy back into the batteries during braking, providing a 10% range extension for greater freedom.

EXPERIENCE INDEPENDENCE.

The innovative design of the twion uses high-tensile aluminum alloys, high-tech plastics and the latest digital motor technology resulting in a sophisticated, light-weight design with easy quick release wheels allowing for smooth, easy transport into a vehicle or public transportation.

twion weighs in at an amazing: 13lbs.*, per wheel including battery!

* Depending on the equipment version.





At a glance

Weight per wheel:	13 lbs.
Speed:	up to 6 mph
Range:	8 - 13 m
Battery technology:	maintenance-free lithium-ion cells
Technologies:	EasyConnex, energy recovery



EXQUISITE DESIGN.

Exclusive, simple and perfectly formed. With every detail, the twion design casts a spell over the beholder. A blend of smooth transitions and distinctive color combinations. High-quality processing and yet intuitive and user-friendly. Alber sets the highest quality standards for the materials used and the design.



ELECTRIFYING.

Another completely new feature is the use of the innovative EasyConnex charging plug. Using magnetic force, it automatically finds the correct plug position. The charging socket on the wheel has an automatic shut off for user safety.

5 LEDs indicate the state of charge of the battery. After only 4 hours of charging, at least 80% of the capacity is again available to you.





REMAIN FLEXIBLE.

twion is suitable for almost every standard active wheelchair. This means that existing, individually adapted wheelchairs can continued to be used as normal. This is made possible by a small, lightweight bracket that is simply fitted to the wheelchair frame.

The user-friendly quick release axle system ensures flexible interchangeability of the manual wheels and the twion wheels. This means that the required drive type is available at all times: manual or with E-Power!



twion[®]
ENERGIZE YOUR WHEELCHAIR



Alber USA LLC
1005 International Drive
Oakdale, PA 15071
Phone (888) 426-8581
Fax (724) 695-3710
info@alber-usa.com
www.alber-usa.com



Scan QR code to receive
further information.