



Invacare® SOLARA™ 3G Adult Tilt-in-Space Wheelchair
AADL Customers Only
AADL # W635
Effective October 1, 2019

NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

FRAME
FRAME TYPE

- SOLARA3G Solara 3G Tilt-in-Space Wheelchair \$3,900
300lb standard weight. A seating system must be purchased for use with this wheelchair. No seating system is provided.

CUSTOM PACKAGES

- HDPG Heavy Duty Package (400 lb. capacity) \$601
Reinforced axles and crossbraces. Weight limit of 400 lbs.

SUSPENSION

- DROCK DeGage Dynamic Rocker Back \$1,456
Not available with HDPG, LOSHEAR, LOSHEARWU, RCANE. The backcane selected will be installed over the DeGage adding approximately 8" to the overall height. See www.degage.us for more details.
- FLSF Frog Legs Suspension Forks \$444
Not available with HDPG.
- SPACK Suspension Package \$1,628
Not available with HDPG, LOSHEAR, LOSHEARWO, RCANE. Includes DROCK and FLSF.

*TRANSPORT TIE-DOWN	
TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lbs to 209 lbs.	TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.
As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.	
* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.	

TRANSPORT TIE DOWN

- Transport Tie-Down options not retrofittable to existing models, not field serviceable.*
- TRBKTS Transport Brackets No Charge
Not available with AVENT, BACKVENT, P750, 2218.

VENT TRAY OPTION

- AVENT Ventilator Tray \$492
Not available with HDPG. Available with 16" wide frame and above and 18" seat-to-floor height and above
- ATBOX Ventilator Battery Box \$171
- BACKVENT Laptop Vent Tray \$1,222
Not available with HDPG. Available with 15" wide frame and above. Not available with LOSHEAR or P750.

TILTING OPTIONS

- THR Traditional Hand Release Standard
- FRLT Cable-Free Foot Release-Left Mount No Charge
- FRRT Cable-Free Foot Release-Right Mount No Charge
- TLL Tilt Lock Out Levers (pair) No Charge
- TLC Tilt Lock Stops No Charge

This is for an additional pair of tilt lock stops, one pair comes standard with chair.

SEAT WIDTH RANGE

- SWR1216 12"-16" Seat Width Range No Charge
Not available with HDPG.
- SWR1418 14"-18" Seat Width Range No Charge
- SWR1620 16"-20" Seat Width Range No Charge
- SWR1822 18"-22" Seat Width Range No Charge
- SWR2024 20" - 24" Seat Width Range No Charge

FRAME WIDTH

- W12 12" Seat Width No Charge
- W13 13" Seat Width No Charge
- W14 14" Seat Width \$159
- W15 15" Seat Width \$159
- W16 16" Seat Width No Charge
- W17 17" Seat Width No Charge
- W18 18" Seat Width No Charge
- W19 19" Seat Width No Charge
- W20 20" Seat Width No Charge
- W21 21" Seat Width \$159
- W22 22" Seat Width \$159
- W23 23" Seat Width \$159
- W24 24" Seat Width \$159

SEAT DEPTH RANGE

- SDR1216 12"-16" Seat Depth Range No Charge
- SDR1519 15"-19" Seat Depth Range No Charge
- SDR1822 18"-22" Seat Depth Range No Charge

Froglegs Suspension Forks

Rear Wheel	Front & Rear Seat-to-Floor Heights				
	Measured from floor to top of seat rail				
	4" 4SPLT	5" 5SOFT 5SPLT 5ALUM 5ALUM2	6" AT902 6ALUM2	6" 1259	8" 1228 AT903
12"	13.50"	15.50"	15.50" 16.00"	16.50" 17.00" 17.50" 18.00"	N/A
16"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00"	N/A	17.00" 17.50" 18.00"	18.50" 19.00"
18"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00"	15.50" 16.00"	18.00"	18.50" 19.00"
20"	15.00"	15.00" 15.50"	15.50" 16.00"	16.50" 17.00" 17.50"	N/A
22"	N/A	N/A	16.00"	16.50" 17.00" 17.50" 18.00"	18.50"
24"	N/A	N/A	N/A	16.50" 17.00" 17.50" 18.00"	18.50" 19.00"

- D12 Seat Depth 12" No Charge
- D13 Seat Depth 13" No Charge
- D14 Seat Depth 14" No Charge
- D15 Seat Depth 15" No Charge
- D16 Seat Depth 16" No Charge
- D15 and D16 are Not Available with FOOTBRK.*
- D17 Seat Depth 17" No Charge
- D18 Seat Depth 18" No Charge
- D19 Seat Depth 19" No Charge
- D20 Seat Depth 20" No Charge
- D21 Seat Depth 21" No Charge
- D22 Seat Depth 22" No Charge

FRAME FINISH

- 24P Wet Black No Charge
- 29P White Snow No Charge
- 30P Sunny Yellow No Charge
- 60P Silver Metallic No Charge
- 61P Electric Red No Charge
- 62P Electric Blue No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge
- 127P Grasshopper Green No Charge

WHEELS

SEAT-TO-FLOOR HEIGHT

- STF125 12.5" Seat-to-Floor No Charge
- STF13 13" Seat-to-Floor No Charge
- STF135 13.5" Seat-to-Floor No Charge
- STF14 14" Seat-to-Floor No Charge
- STF145 14.5" Seat-to-Floor No Charge
- STF15 15" Seat-to-Floor No Charge
- STF155 15.5" Seat-to-Floor No Charge
- STF16 16" Seat-to-Floor No Charge
- STF165 16.5" Seat-to-Floor No Charge
- STF17 17" Seat-to-Floor No Charge
- STF175 17.5" Seat-to-Floor No Charge
- STF18 18" Seat-to-Floor No Charge
- STF185 18.5" Seat-to-Floor No Charge
- STF19 19" Seat-to-Floor No Charge

CASTER OPTIONS

- 4SPLT 4 Composite Urethane Light Up No Charge
- Not available with HDPG.*
- 4ALUM2 4" x 1.4 Aluminum Soft Roll \$134
- 51DC 5" x 1" Urethane Casters No Charge
- Not available with HDPG.*
- 5SPLT 5" x 1" Composite Urethane Light Up No Charge
- Not available with HDPG.*
- 5ALUM2 5" x 1.4 Aluminum Soft Roll \$134
- 61DC 6 in x 1 in Endura Roll Caster (Urethane) No Charge
- Not available with HDPG.*
- 1259 6" X 2" Semi-Pneumatic Caster \$132
- Not available with HDPG.*
- 815DC 8 in x 1.5 in Endura Roll Caster (Urethane) \$150
- 6ALUM2 6" x 1.4 Aluminum Soft Roll \$134
- 81DC 8 in x 1 in Endura Roll Caster (Urethane) No Charge

WHEEL OPTIONS

- 12COM 12" Composite **Standard**
- 16COM 16" Composite Wheels No Charge
- Not available with urethane tires. Only available with aluminum handrims.*
- 18COM 18" Composite Wheels No Charge
- Not available with urethane tires. Only available with aluminum handrims.*
- 20COM 20" Composite No Charge
- 22COM 22" Composite No Charge
- 24COM7 24" Composite-7 Spoke No Charge

AXLE TYPE

- A1577 Quick Release Axles No Charge
- PAXLE Permanent Axles No Charge

TIRE OPTIONS

- PNEU Pneumatic No Charge
- Not available with HDPG.*
- PNEUFF Pneumatic-Flat Free \$128
- URETH Urethane (Treaded) No Charge

Front and Rear Seat-to-Floor Heights (measured from Floor to Top of Seat Rail)					
Rear Wheel	4" 4SPLT 4ALUM2	5" 51DC 5ALUM2 5SPLT	6" 61DC 6ALUM2	6" 1259	8" 81DC 815DC
12" Pneu PneuFF	12 1/2" 13" 13 1/2" 15 1/2" 16"	13 1/2" 15 1/2" 16" 16 1/2"	15 1/2" 16" 16 1/2" 17"	15 1/2" 16" 16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18"
12" Urethane	12 1/2" 14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	16 1/2" 17"
16"	12 1/2" 13" 13 1/2" 14" 14 1/2" 15"	13 1/2" 14" 14 1/2" 15"	14 1/2" 15" 17"	14 1/2" 15" 17 1/2" 18"	17" 17 1/2" 18 1/2" 19"
18"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16" 16 1/2" 18"	18" 18 1/2" 19"
20"	15" 15 1/2" 16"	15" 15 1/2" 16" 16 1/2"	15 1/2" 16" 16 1/2" 17"	15 1/2" 16" 16 1/2" 17" 17 1/2"	16 1/2" 17" 17 1/2"
22"	16"	16" 16 1/2"	16" 16 1/2" 17"	16" 16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18" 18 1/2"
24"	N/A	16 1/2"	16 1/2" 17"	16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18" 18 1/2" 19"

All Dimensions Are Approximate

HANDRIM OPTIONS

- AHDRM Aluminum Handrims No Charge
- CHDRM Chrome Plated Handrims No Charge
- 1512 Plastic Coated Aluminum Handrims \$114
- AT1224-8 Plastic Coated w/8 Oblique Projections \$307
Only available with 24COM7.
- AT1225-8 Plastic Coated w/8 Vertical Projections \$307
Only available with 24COM7.
- AT1225-12 Plastic Coated w/12 Vertical Projections \$307
Only available with 24COM7.
- L8702 Omit Handrims Less net(\$15)

Seat Pan Growth Range (1" increments)	
Chair Width	Seat Pan Width Range
W12	12"-14"
W13	12"-14"
W14	13"-16"
W15	13"-16"
W16	15"-18"
W17	15"-18"
W18	17"-20"
W19	17"-20"
W20	19"-22"
W21	19"-22"
W22	21"-24"
W23	21"-24"
W24	23"-24"

WHEEL LOCKS

- FOOTBRK Foot Operated Push To Lock Wheel Lock No Charge
Not available with AVENT(vent tray)
- PULLOCK Pull To Lock Wheel Locks No Charge
- PUSHLOCK Push To Lock Wheel Locks No Charge
- EXTLOCK Push To Lock Wheel Locks w/6" Bolt-On Extensions \$64
- FOOTHUB Foot Operated Hublock \$268
Not available with 16COM
- ATHUB Cane Mount Hub Lock No Charge
Not available with 16COM. Not available with LOSHEAR.
- USEHUB Seat Rail Mount Hub Lock No Charge
Not available with 16COM. Only available with TAD or TAF.

POSITIONING

SEAT STYLE

Both pans grow in width and require no new hardware or parts when converting from flush to drop or vice versa. They do not come with the chair and thus, must be selected if needed.

- FLSEAT Flush Mounted Seat Pan \$132
- DPSEAT Drop Mounted Seat Pan \$132

ARM HEIGHT

- SARMHT Low Arm Height No Charge
- HARMHT High Arm Height No Charge
Not available with LOSHEAR, LOSHEARWU or Cantilever arms.

Arm Heights (measured from Top of Seat Rail to Top of Armpad)		
Arm Type	Low	High
TAD	6 1/2" - 12 1/2"	9" - 17"
TAF	1/2" increments	1" increments
FXDESK	9 1/2"	12"
FXFULL		
ADJDESK	9 1/2" - 13 1/2"	12" - 16"
ADJFULL	1" increments	1" increments
FADJDESK	9" - 13"	12" - 16"
FADJFULL	1" increments	1" increments

ARM TYPE (PAIRS)

Two Point arms have been re-designed.

- FXDESK Desk Length Fixed Height Conv \$261
Available on 15" seat depth and above.
- FXFULL Full Length Fixed Height Conv \$261
Available on 15" seat depth and above.
- ADJDESK Desk Length Adjustable Height \$372
Available on 15" seat depth and above.
- ADJFULL Full Length Adjustable Height \$372
Available on 15" seat depth and above.
- FADJDESK Two-Point Flip-Back Adjustable Ht Desk \$372
Not available with HDPG. Available on 15" seat depth and above. Not available with LOSHEAR, LOSHEARWO.
- FADJFULL Two-Point Flip-Back Adjustable Ht Full \$372
Not available with HDPG. Available on 15" seat depth and above. Not available with LOSHEAR, LOSHEARWO.
- TAF Full Adjustable Height \$372
- TAD Desk Adjustable Height \$372

Back Cane Heights (measured from Top of Seat Rail)		
Back Cane Type	Inches to 10 degree bend	Total Height
BH20	Straight Cane	20.5"
BH20B/BH24B	15 1/2"	24.5"
BH24	Straight Cane	24 1/2"
RECLINER	Straight Cane	26"
RECLINER + Headrest	Straight Cane	36"
ABH1720	16.2" - 21.2"	29.3" - 34.3"

CANTILEVER ARMS

Not available with LOSHEAR, LOSHEARWU OR TRBKTS.

- AT99N Cantilever Non Locking Bolt On Arms \$261
- AT55D Cantilever Locking Desk Lgth.Bolt On Arms \$372
- AT55F Cantilever Locking Full Lgth.Bolt On Arms \$372

BACK TYPE

- SCANE Standard Aluminum Back Cane **Standard**
- RCANE Steel-Reinforced Back Cane No Charge
- LOSHEAR Low Shear Reclining Back with Upholstery \$1,232
Not Available on HDPG. Comes with back upholstery and AASH. Available on 14" - 20" seat width.
- LOSHEARWO Low Shear Reclining Back without Upholstery \$1,212
Not available with HDPG. Not available with a headrest and pillow. Comes standard with AASH.
- RECEXT Reclining Headrest Pillow and Extension \$131
- RLL Recliner Lockout Levers No Charge
Not available with SCANE or RCANE.
- RLC Recliner Lock Stops \$31

BACK HEIGHT

- BH20 Back Height 20" No Charge
Not available with P750.
- BH20B Back Height 20"with bend No Charge
Not available with P750.
- BH24 Back Height 24" **Standard**
- BH24B Back Height 24"with bend No Charge
- ABH1720 Adjustable Height Back Canes (17" - 20") No Charge

Not available with P750. Must choose the AASH when ordering the ABH1720. The ABH1720 is not available with the 1483 Non-Adjustable Rigidizing Handle Bar. Also not available with RCANE because it already consists of heavy duty construction.

FOLD OPTION BACK

<input type="checkbox"/> FXBACK	Non-Fold Down Back Canes	Standard
<input type="checkbox"/> FLDBACK	Fold-Down Back Canes	No Charge

Not available with LOSHEAR

SEAT STRAPS

<input type="checkbox"/> 1311BK	Airline Buckle Seat Positioning Strap	\$84
<input type="checkbox"/> 1515	Auto Style Buckle Seat Positioning Strap	\$84
<input type="checkbox"/> 1321BK	Hook and Loop Seat Positioning Strap	\$66

CHEST STRAPS

<input type="checkbox"/> 7311BK	Airline Buckle Chest Positioning Strap	\$84
<input type="checkbox"/> 9515	Auto Style Chest Positioning Strap	\$84
<input type="checkbox"/> 7321BK	Hook and Loop Chest Positioning Strap	\$66

BODYPOINT® PADDED CALF STRAPS

See www.bodypoint.com for complete details. Products are not installed on chair. Measure chair width to determine size

<input type="checkbox"/> BPCALFS	Padded Calf Straps (13"-15")	\$54
<input type="checkbox"/> BPCALFM	Padded Calf Straps (15"-18")	\$58
<input type="checkbox"/> BPCALFL	Padded Calf Straps (18"-22")	\$67

BODYPOINT® ANKLE HUGGERS®

See www.bodypoint.com for complete details. Products are not installed on chair. Measure circumference above ankle to determine size.

<input type="checkbox"/> BPAHGS	Side-Release Ankle Huggers (6.5"-8")	\$162
<input type="checkbox"/> BPAHGM	Side-Release Ankle Huggers (7.5"-9")	\$166
<input type="checkbox"/> BPAHGL	Side-Release Ankle Huggers (8.5"-11.5")	\$166
<input type="checkbox"/> BPAHGXL	Side-Release Ankle Huggers (11"-13")	\$171

BODYPOINT® STAYFLEX™

See www.bodypoint.com for complete details. Products are not installed on chair. Measure shoulder width to determine size.

BPSFM and BPSFL are approx 2" wider in center than BPSFNM and BPSFNL.

<input type="checkbox"/> BPSFM	Anterior Trunk Support w/Zip (13"-16")	\$181
<input type="checkbox"/> BPSFL	Anterior Trunk Support w/Zip (16"-19")	\$185
<input type="checkbox"/> BPSFNM	Anterior Trunk Support Narrow w/Zip (13"-16")	\$181
<input type="checkbox"/> BPSFNL	Anterior Trunk Support Narrow w/Zip (16"-19")	\$185

MISCELLANEOUS

HEADREST ACCESSORIES (OFF CHAIR)

Headrests will ship separately from chair. Non standard chair discounts apply. Please call for more information.

<input type="checkbox"/> 5500071	Comfort Headrest Assembly w/ Multi Axis Hardware	\$375
<input type="checkbox"/> 5500061	Basic Headrest Kit	\$279
<input type="checkbox"/> 1130083	Comfort Headrest w/ Standard Mounting Hardware	\$198

BACK ACCESSORIES

Only available with TAD or TAF. Not available with LOSHEAR.

<input type="checkbox"/> ADJRGD	Adjustable Width Back Rigidizing Bar	\$48
<input type="checkbox"/> 1483	Non-Adjustable Stroller Handles	\$156
<input type="checkbox"/> AASH	Adj. Angle Stroller Handles (Outside Mount)	\$454
<input type="checkbox"/> AASHIN	Adj. Angle Stroller Handles (Inside Mount)	\$454

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> A1580	Rear Anti-Tippers	Standard
<input type="checkbox"/> BP02	Backpack	No Charge
<input type="checkbox"/> 1209	Tire Pump	\$63
<input type="checkbox"/> P750	Oxygen Holder	\$163

Not available with HDPG. Not available with LOSHEAR, LOSHEARWO OR TRBKTS. Not available with BH20 or ABH1720.

Not available with AVENT (vent tray). Not available with 8" casters. Not available on seat widths less than 16".

Not available with LOSHEAR, LOSHEARWO OR TRBKTS.

FRONT RIGGING OPTIONS


FRONT RIGGING MOUNTING

<input type="checkbox"/> STD MNT	STD Front Rigging Mount	Standard
<input type="checkbox"/> OFF MNT	Offset Front Rigging Mount	No Charge

OFF MNT re-designed for a smoother and better angle. Not available with 12"-16" seat depth range

FRONT RIGGING

<input type="checkbox"/> 70LIFTNT	70° Lift Off Non-Taper Footrests	\$241
<input type="checkbox"/> 70LIFTMFX	70° Lift Off MFX Footrests	\$376
<i>Must select AT5543</i>		
<input type="checkbox"/> 90LIFTMFX	90° Lift Off MFX Footrests	\$248
<i>Must select AT5543</i>		
<input type="checkbox"/> SWING1319	60° Swingaway Footrest	\$241
<input type="checkbox"/> SWING914	60° Swingaway Footrest MFX	\$376
<input type="checkbox"/> 70NTAPER	70° Swingaway Footrest	\$241
<input type="checkbox"/> 70MFX	70° Swingaway Footrest MFX	\$376
<i>Must select AT5543</i>		
<input type="checkbox"/> 70STAPER	70° Tapered Swingaway Footrest	\$281
<i>Available on seat widths 16" and above</i>		
<input type="checkbox"/> 80NTAPER	80° Non-Tapered Swingaway Footrest	\$281
<input type="checkbox"/> 80MFX	80° Swingaway Footrest MFX	\$376
<i>Must select AT5543</i>		
<input type="checkbox"/> 90SWING	90° Swingaway Footrest	\$248
<i>Must select AT5543</i>		
<input type="checkbox"/> AT5044LIFT	Elevating Legrest	\$567
<input type="checkbox"/> AT5044	Elevating Legrest	\$567
<input type="checkbox"/> BILAT	Bilateral Contracture Footrest	\$457
<i>Not available with HDPG.</i>		
<input type="checkbox"/> ATL4	Centre Pivot Articulating Legrests	\$392



Available Through Our Specials Departments

- 20" High Reclining Back SM-050503023551
- 25" - 26" Wide SM-010312110036
- 80° Footrest Hanger SM-200411121756
- Front Rigging Mount: Standard on one side and offset on the opposite side C-090731110749

Not available with 1651.

110CP Contracture Platform

Not available with HDPG. Available on 16"- 22" seat widths.

\$457



FOOTPLATE TYPE

- 1651 Flip-Up Composite Footplate No Charge
Not available with Bodypoint Ankle Huggers.
- A1562 Aluminum Footplate No Charge
Not available with Bodypoint Ankle Huggers.
- A1564 Extra Large Aluminum Footplate \$106
Not available with Bodypoint Ankle Huggers.
- AT5543 Adjustable Angle Flip-Up Footplates \$132

ADJUSTABLE ANGLE FOOTPLATES RIGHT

Must select with BILAT.

(FOR BI-LATERAL CONTRACTURAL FOOTRESTS ONLY)

- AAF325R Adj. Angle Footplate Right 3.25" No Charge
- AAF425R Adj. Angle Footplate Right 4.25" No Charge
- AAF525R Adj. Angle Footplate Right 5.25" No Charge
- AAF625R Adj. Angle Footplate Right 6.25" No Charge
- AAF725R Adj. Angle Footplate Right 7.25" No Charge

ADJUSTABLE ANGLE FOOTPLATES LEFT

Must select with BILAT.

(FOR BI-LATERAL CONTRACTURAL FOOTRESTS ONLY)

- AAF325L Adj. Angle Footplate Left 3.25" No Charge
- AAF425L Adj. Angle Footplate Left 4.25" No Charge
- AAF525L Adj. Angle Footplate Left 5.25" No Charge
- AAF625L Adj. Angle Footplate Left 6.25" No Charge
- AAF725L Adj. Angle Footplate Left 7.25" No Charge

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$63
- CLFPDLRG Large Calf Pad \$268
Only available on Elevating leg rests.
- 9342 H-Heel Strap \$94
- FRMPRBK ChairArmor Black Neoprene Frame Protectors No Charge
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- FRMPRR ChairArmor Red Neoprene Frame Protectors \$70
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- FRMPRB ChairArmor Blue Neoprene Frame Protectors \$70
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- ALPT Longer Pivot And Slide Tube \$132
- ALPT3 3" Extension 90° footrests \$132
- SCLAMP Clamp-On Sector Blocks \$57
Not available with 70LIFTNT, 70LIFTMFX, 90LIFTMFX, AT5044LIFT, BILAT, 110CP. Adjustable angle. See Front Rigging Mounts image.

HEEL LOOP/TOE LOOP

- HLOOP Heel Loops Without Ankle Straps \$79
- 1345 Toe Loops \$64

Knee-to-Heel Range (Measured from edge of seat pan to back of footplate)		
Front Rigging Type	Minimum	Maximum
Range adjusts in 1/2" increments.		
60°	13 1/2"	17 1/2"
60° MFX	10"	13 1/2"
*70°	13 1/2"	18 1/2"
*70° Tapered	13 1/4"	17.00"
80°	10.5"	16"
Bilateral Contracture STF125-STF16	5 1/2"	9"
Bilateral Contracture STF165-STF19	5 1/2"	13"
Range adjusts in 1" increments.		
70° MFX	6 1/2"	12 1/2"
80° MFX	4.5"	9"
90°	4 1/2"	7 1/2"
Infinite Range Adjustments		
*Elevating	14.00"	20.00"
Articulating	16.00"	24.00"
Contracture Platform STF125-STF16	7 1/4"	8 1/4"
Contracture Platform STF165-STF19	11 1/2"	16"
* ALPT shifts range by 2"		

Specifications	
Frame Style	Tilt-in-Space
Overall Width	Seat width + 8 1/4" (12" wheels) Seat width + 9 1/4" (18-24" wheels)
Seat Widths	12" - 24"
Seat Depths	12" - 22"
Back Height	17" - 20", 20", 24"
Adjustable Angle Back	0°, 5°, 10°, 15°, 20°, 25°, 30°
Tilt Range	-5° to 50°
Recliner Back	175°
Weight Limit	300 lbs./400 lbs. with Heavy-Duty Package
Weight	34 1/2 lbs*
Warranty	Lifetime on frame and crossbraces
*Base chair without Accessories	

Tilt Assist Explained
The Solara® 3G wheelchair is shipped from the factory pre-configured with one of three "Tilt Assist" settings: *Front, Mid or Rear*, based solely on the seat depth range and seat depth ordered. This setting describes how the user's weight is distributed on the frame itself, i.e., does one front load the chair, rear load or is the weight evenly distributed. As compared to the 2G, the 3G requires up to 30% less force needed to tilt because of these settings. If factors other than seat depth or seat depth range come into play that affect the weight distribution of the chair, the provider has the ability to change the Tilt Assist setting in the field by simply adjusting 2 screws on each side of the frame.

Tilt Assist Factory Setting		
	Seat Depth Range 15" - 19"	Seat Depth Range 18" - 22"
D15	Front Setting	
D16	Front Setting	
D17	Mid Setting	
D18	Rear Setting	Front Setting
D19	Rear Setting	Front Setting
D20		Mid Setting
D21		Rear Setting
D22		Rear Setting


The factory setting may need to be adjusted if a seating system is added. The factory setting does not account for conditions outside of the seat depth and seat depth range which may cause more weight to be distributed to the front or back of the chair.

When adjusting the Tilt Assist there are 4 bolts to change (2 on each side.) See the illustration to the left and chart above.

The **Rear** setting has the rocker tube telescoped entirely into the seat frame support tube where the adjustment bolt is in the bottom hole. Along the base frame the bolt is in the top hole. (e.g. accommodates a seat depth that is in the highest end of the seat depth range)

For the **Mid** setting the adjustment bolt is simply in the middle hole along the rocker tube which telescopes out and the middle hole along the base frame. (e.g. optimal for a seat depth in the middle of the seat depth range)

The **Front** setting has the bolt in the top hole of the rocker tube which is fully telescoped out of the seat frame support tube and in the bottom hole along the base frame. (e.g. accommodates a seat depth in the lowest end of the seat depth range)



Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)
Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

WWW.INVACARE.CA

NOTE: All specifications and dimensions are approximate
For the most current pricing information visit www.invacare.ca