



# Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

## Thank you for using Invacare eForms.

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit [www.invacare.com/eForms](http://www.invacare.com/eForms)



## Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

### Save

---

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

### Submit

---

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

### Print

---

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**800-678-4682**



## Account Information

Request Type:      Quote                      Order

Date: \_\_\_\_\_

Account #: \_\_\_\_\_

Company: \_\_\_\_\_

SHIP TO  
Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Country: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

CONTACT  
Name: \_\_\_\_\_

Back Up Contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Comments: \_\_\_\_\_

## Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

## FRAME SELECTIONS

### FRAME TYPE

SOLARA3G Solara 3G Tilt-in-Space wheelchair .....\$2,982.00  
E1161

A seating system must be purchased. No seating system is provided.  
300 lb standard weight capacity. See www.invacare.com/Solara3G  
for additional information when configuring this chair.

### CUSTOM PACKAGES

HDPG Heavy Duty Package (400 lbs. Weight Capacity) .....  
.....\$400.00 K0108

Reinforced axles and crossbraces. Weight limit of 400 lbs.

## SUSPENSION SELECTION

### SUSPENSION

DROCK Seating Dynamics Rocker Back.....\$1985.00 E2398

Not available with HDPG, LOSHEAR, LOSHEARWO, RCANE, FLDBACK.  
See www.seatingdynamics.com for more details.

FLSF Frog Legs Suspension Forks.....\$398.00 E1015

Not available with HDPG. See chart page 7.

SPACK Suspension Package .....\$1449.00 E1015

Not available with HDPG, LOSHEAR, LOSHEARWO, RCANE, FLDBACK.  
Includes DROCK and FLSF.

### Transport Tie Down

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tiedown systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

## TRANSPORT/VENT TRAY SELECTION

### TRANSPORT TIE DOWN

TRBKTS Transport Brackets .....\$300.00 K0108

Not available with AVENT, BACKVENT, P750, 2218. Now available as a service part kit for existing models.

### VENT TRAY OPTION

AVENT Ventilator Tray (each).....\$295.00 E1029

Not available with HDPG. Available with 16" wide frame and above and 18" Seat-to-Floor Height and above.

ATBOX Ventilator Battery Box (each) .....\$250.00

BACKVENT Laptop Vent Tray .....\$1,540.00 E1030

Not available with HDPG, LOSHEAR, P750. Available with 15" wide frame and above.



# FRAME SELECTIONS

## TILTING OPTIONS

- FRLT Cable-Free Foot Release-Left Mount..... No Charge
- FRRT Cable-Free Foot Release-Right Mount.....No Charge
- THR Traditional Hand Release ..... Standard

- TLL Tilt Lock Levers..... No Charge
- TLC Tilt Lock Stops (pair)..... \$35.00  
Additional pair of tilt lock stops, one pair comes standard with chair.

## SEAT WIDTH RANGE

- SWR1216 12" – 16" Seat Width Range ..... No Charge  
Not available with HDPG.
- SWR1418 14" – 18" Seat Width Range..... No Charge
- SWR1620 16" – 20" Seat Width Range ..... No Charge

- SWR1822 18" – 22" Seat Width Range ..... No Charge
- SWR2024 20" – 24" Seat Width Range..... No Charge

## SEAT WIDTH

- W12 12" Wide ..... No Charge
- W13 13" Wide ..... No Charge
- W14 14" Wide ..... No Charge
- W15 15" Wide ..... No Charge
- W16 16" Wide ..... No Charge
- W17 17" Wide ..... No Charge
- W18 18" Wide ..... No Charge

- W19 19" Wide ..... No Charge
- W20 20" Wide ..... \$500.00 E2201
- W21 21" Wide ..... \$500.00 E2201 K0108
- W22 22" Wide ..... \$500.00 E2201 K0108
- W23 23" Wide ..... \$500.00 E2201 K0108
- W24 24" Wide ..... \$500.00 E2202 K0108

## SEAT DEPTH RANGE

- SDR1216 12" – 16" Seat Depth Range..... No Charge  
Tilt Assist feature not available on this frame.
- SDR1519 15" – 19" Seat Depth Range ..... No Charge

- SDR1822 18" – 22" Seat Depth Range..... No Charge

## SEAT DEPTH

- D12 12" Seat Depth ..... No Charge
- D13 13" Seat Depth ..... No Charge
- D14 14" Seat Depth ..... No Charge
- D15 15" Seat Depth ..... No Charge
- D16 16" Seat Depth ..... No Charge
- D17 17" Seat Depth ..... No Charge

- D18 18" Seat Depth ..... No Charge
- D19 19" Seat Depth ..... No Charge
- D20 20" Seat Depth..... \$550.00 E2203
- D21 21" Seat Depth..... \$550.00 E2203 K0108
- D22 22" Seat Depth..... \$900.00 E2204 K0108



## FRAME FINISH

### FRAME COLORS

- 24P Wet Black ..... No Charge
- 29P White Snow ..... No Charge
- 60P Silver Metallic ..... No Charge
- 115P Black Prism ..... No Charge
- 105P Midnight Blue ..... No Charge
- 116P Sunkissed Orange ..... No Charge
- 117P Tangerine Red ..... No Charge
- 30P Sunny Yellow ..... No Charge

- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge
- 62P Electric Blue ..... No Charge
- 127P Grasshopper Green ..... No Charge
- Touch Up Paint ..... \$20.00

Touch up paint is not available in all colors.

### Touch Up Color

Discontinued frame colors are available through our specials department. C-070109030544

## CASTER, WHEEL AND TIRE SELECTIONS

### SEAT-TO-FLOOR HEIGHT

- STF125 Seat-To-Floor 12 1/2" (10 1/2" w/DPSEAT)..... No Charge
- STF13 Seat-To-Floor 13" (11" w/DPSEAT)..... No Charge
- STF135 Seat-To-Floor 13 1/2" (11 1/2" w/DPSEAT) ..... No Charge
- STF14 Seat-To-Floor 14" (12" w/DPSEAT) ..... No Charge
- STF145 Seat-To-Floor 14 1/2" (12 1/2" w/DPSEAT)..... No Charge
- STF15 Seat-To-Floor 15" (13" w/DPSEAT) ..... No Charge
- STF155 Seat-To-Floor 15 1/2" (13 1/2" w/DPSEAT)..... No Charge

- STF16 Seat-To-Floor 16" (14" w/DPSEAT) ..... No Charge
- STF165 Seat-To-Floor 16 1/2" (14 1/2" w/DPSEAT)..... No Charge
- STF17 Seat-To-Floor 17" (15" w/DPSEAT) ..... No Charge
- STF175 Seat-To-Floor 17 1/2" (15 1/2" w/DPSEAT)..... No Charge
- STF18 Seat-To-Floor 18" (16" w/DPSEAT) ..... No Charge
- STF185 Seat-To-Floor 18 1/2" (16 1/2" w/DPSEAT) ..... No Charge
- STF19 Seat-To-Floor 19" (17" w/DPSEAT) ..... No Charge

### CASTER OPTIONS (PAIR)

- 4SPLT 4" x 1" Composite Urethane Light Up..... No Charge  
Not available with HDPG.
- 4ALUM2 4" x 1.4" Aluminum Soft Roll ..... \$125.00
- 5SPLT 5" x 1" Composite Urethane Light Up..... No Charge  
Not available with HDPG.
- 51DC 5" x 1" Urethane Casters..... No Charge  
Not available with HDPG.
- 5ALUM2 5" x 1.4" Aluminum Soft Roll ..... \$125.00

- 61DC 6" x 1" Urethane Caster..... No Charge  
Not available with HDPG.
- 1259 6" x 2" Semi-Pneumatic Casters..... \$90.00  
Not available with HDPG.
- 6ALUM2 6" x 1.4" Aluminum Soft Roll ..... \$125.00
- 81DC 8" x 1" Urethane Casters ..... No Charge
- 815DC 8" x 1.5" Urethane Caster ..... No Charge



# Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

## WHEEL OPTIONS (PAIR)

- 12COM 12" Composite Wheels ..... No Charge
- 16COM 16" Composite Wheels ..... No Charge  
Not available with urethane tires. Only available with aluminum handrims.
- 18COM 18" Composite Wheels ..... No Charge  
Not available with urethane tires. Only available with aluminum handrims.

- 20COM 20" Composite Wheels ..... No Charge
- 22COM 22" Composite Wheels ..... No Charge
- 24COM 24" Composite Wheels ..... No Charge

## AXLE TYPE (PAIR)

- A1577 Quick Release Axles ..... No Charge

- PAXLE Permanent Axles ..... Standard

## TIRE OPTIONS (PAIR)

- PNEU Pneumatic Tires ..... \$85.00  
Not available with HDPG.
- PNEUFF Flat Free Inserts ..... \$85.00 E2213
- URETH Urethane (Treaded) Tires ..... No Charge

## HANDRIM OPTIONS (PAIR)

- CHDRM Chrome Plated Handrims ..... No Charge
- 1512 Plastic Coated Aluminum Handrims ..... \$115.00
- AT-1225-8 Plastic Coated with 8 Vertical Projections ..... \$295.00  
Only available with 24COM7.
- AT-1224-8 Plastic Coated with 8 Oblique Projections ..... \$295.00  
Only available with 24COM7.
- AT1225-12 Plastic Coated with 12 Vertical Projections ..... \$295.00  
Only available with 24COM7.
- L8702 Omit Handrims ..... \$0.00

Front and Rear Seat-to-Floor Heights (measured from Floor to Top of Seat Rail)					
Rear Wheel	4" 4SPLT 4ALUM2	5" 51DC 5ALUM2 5SPLT	6" 61DC 6ALUM2	6" 1259	8" 81DC 815DC
12" Pneu PneuFF	12 1/2" 13" 13 1/2" 15 1/2" 16"	13 1/2" 15 1/2" 16" 16 1/2"	15 1/2" 16" 16 1/2" 17"	15 1/2" 16" 16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18"
12" Urethane	12 1/2" 14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	16 1/2" 17"
16"	12 1/2" 13" 13 1/2" 14" 14 1/2" 15"	13 1/2" 14" 14 1/2" 15"	14 1/2" 15" 17"	14 1/2" 15" 17" 17 1/2" 18"	17" 17 1/2" 18" 18 1/2" 19"
18"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16" 18"	18" 18 1/2" 19"
20"	15" 15 1/2" 16"	15" 15 1/2" 16" 16 1/2"	15" 15 1/2" 16" 16 1/2" 17"	15" 15 1/2" 16" 16 1/2" 17" 17 1/2"	16 1/2" 17" 17 1/2"
22"	16"	16" 16 1/2"	16" 16 1/2" 17"	16" 16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18" 18 1/2"
24"	N/A	16 1/2"	16 1/2" 17"	16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18" 18 1/2" 19"

All Dimensions Are Approximate



## WHEEL LOCKS

- FOOTBRK Foot Operated Push To Lock Wheel Lock ..... \$75.00 K0108  
Not available with AVENT.
- PULLLOCK Pull To Lock Wheel Locks (pair) ..... No Charge
- PUSHLOCK Push To Lock Wheel Locks (pair).....No Charge
- EXTLOCK Push To Lock with 6 Bolt-On Ext (pair).....\$75.00  
E0961

- FOOTHUB Foot Operated Hublock ..... \$200.00 K0108  
Not available with 16COM.
- ATHUB Cane Mount Hub Lock ..... \$200.00 K0108  
Not available with LOSHEAR or LOSHEARWO, 16COM.
- USEHUB Seat Rail Mount Hub Lock ..... \$200.00 K0108  
Only available with TAD of TAF. Not available with 16COM.

## ARM AND BACK SELECTIONS

### SEAT STYLE

FLSEAT or DPSEAT are solid seat pans and do not come standard with chair. Both pans grow in width and require no new hardware or parts when converting from flush to drop or vice versa.

- FLSEAT Flush Seat Pan ..... \$160.00 E2231
- DPSEAT Drop Seat Pan ..... \$160.00 E2231  
Accommodates 2" lower seat to floor heights. Use a cushion not less than 3" thick.

Arm Heights (measured from Top of Seat Rail to Top of Armpad)		
Arm Type	Low	High
TAD TAF	6 1/2" - 12 1/2" 1/2" increments	9" - 17" 1" increments
FXDESK FXFULL	9 1/2"	12"
ADJDESK ADJFULL	9 1/2" - 13 1/2" 1" increments	12" - 16" 1" increments
FADJDESK FADJFULL	9" - 13" 1" increments	12" - 16" 1" increments

Seat Pan Growth Range (1" increments)	
Chair Width	Seat Pan Width Range
W12	12"-14"
W13	12"-14"
W14	13"-16"
W15	13"-16"
W16	15"-18"
W17	15"-18"
W18	17"-20"
W19	17"-20"
W20	19"-22"
W21	19"-22"
W22	21"-24"
W23	21"-24"
W24	23"-24"

### ARM HEIGHT

- SARMHT Low Arm Height..... No Charge

- HARMHT High Arm Height..... No Charge

### CANTILEVER ARMS

Not available with LOSHEAR or LOSHEARWO

- AT99N Cantilever Non Locking Bolt on Arms .....No Charge
- AT55D Cantilever Locking Desk Length Bolt on Arms ..... \$265.00

- AT55F Cantilever Locking Full Length Bolt on Arms..... \$265.00

### T-ARMS (PAIR)

- TAD Adjustable Height, Desk Length .....\$285.00 E0973

- TAF Adjustable Height, Full Length .....\$285.00 E0973

### TWO-POINT, FIXED HEIGHT ARMS (PAIR)

Available on 15" seat depth and above.

- FXDESK Conventional, Desk Length .....No Charge

- FXFULL Conventional, Full Length .....No Charge



**TWO-POINT, ADJUSTABLE HEIGHT ARMS (PAIR)**

Available on 15" seat depth and above.

ADJDESK Conventional, Desk Length ..... \$265.00 E0973

ADJFULL Conventional, Full Length ..... \$265.00 E0973

FADJDESK Flip-Back, Desk Length ..... \$295.00 E0973

Not available with LOSHEAR, LOSHEARWO or HDPG.

FADJFULL Flip-Back, Full Length ..... \$295.00 E0973

Not available with LOSHEAR, LOSHEARWO or HDPG.

**BACK TYPE**

SCANE Standard Aluminum Back Cane ..... Standard

RCANE Steel-Reinforced Back Cane ..... No Charge

LOSHEAR Low Shear Reclining Back with Upholstery .....  
..... \$910.00 E1226

Not available with HDPG. Comes with back upholstery and AASH. Available on 14"-20" seat width.  
Not available with Foot Operated Wheel Lock

LOSHEARWO Low Shear Reclining Back without Upholstery .....  
..... \$880.00 E1226

Not available with HDPG, headrest and pillow. Comes standard with AASH.  
Not available with Foot Operated Wheel Lock

RECEXT Reclining Headrest Pillow and Extension .....  
..... \$85.00 E0966

RLL Recliner Lockout Levers ..... No Charge

Not available with SCANE or RCANE.

RLC Recliner Lock Stops ..... \$35.00

Back Cane Heights (measured from Top of Seat Rail)		
Back Cane Type	Inches to 10 degree bend	Total Height
BH20	Straight Cane	20.5"
BH20B/BH24B	15 1/2"	24.5"
BH24	Straight Cane	24 1/2"
LOSHEAR	Straight Cane	26"
LOSHEAR + Headrest	Straight Cane	36"
ABH1720	16.2" - 21.2"	29.3" - 34.3"

**BACK HEIGHT**

BH20B Back Height 20" with bend ..... No Charge  
Not available with P750.

BH20 Strt Cane w Push Handle Back Ht 20" ..... No Charge  
A9900  
Not available with P750.

BH24B Back Height 24" with bend ..... No Charge

BH24 Back Height 24" ..... No Charge

ABH1720 Adj Height Back Canes (17" – 20") ..... No Charge  
Not available with RCANE (because it consists of heavy duty construction), 1483, ATHUB, P750.

**BACK CANE OPTIONS**

FXBACK Non-Fold Down Back Canes ..... Standard

FLDBACK Fold-Down Back Canes ..... No Charge  
Not available with LOSHEAR, DROCK.

**POSITIONING STRAPS (EACH)**

1515 Push Button Buckle Seat Positioning Strap ..... \$60.00  
E0978

1321BK Hook and Loop Seat Positioning Strap ..... \$60.00  
E0978

9515 Push Button Chest Positioning Strap ..... \$117.00 E0960  
Not available with LOSHEAR or LOSHEARWO.

7321BK Hook and Loop Chest Positioning Strap ..... \$117.00  
E0960  
Not available with LOSHEAR or LOSHEARWO.



### BODYPOINT® PADDED CALF STRAPS

See [www.bodypoint.com](http://www.bodypoint.com) for complete details. Measure chair width to determine size.

BPCALFS Padded Calf Straps (13" - 15").....	\$68.00 K0038	BPCALFL Padded Calf Straps (18" - 22") .....	\$68.00 K0038
BPCALFM Padded Calf Straps (15" - 18") .....	\$68.00 K0038		

### BODYPOINT® ANKLE HUGGERS®

See [www.bodypoint.com](http://www.bodypoint.com) for complete details. Measure circumference above ankle to determine size.

BPAHGS Side - Release Ankle Huggers (6.5" - 8").....	\$180.00 K0108	BPAHGL Side - Release Ankle Huggers (8.5" - 11.5") .....	\$180.00 K0108
BPAHGM Side - Release Ankle Huggers (7.5" - 9").....	\$180.00 K0108	BPAHGXL Side - Release Ankle Huggers (11" - 13") .....	\$180.00 K0108

### BODYPOINT® STAYFLEX™

See [www.bodypoint.com](http://www.bodypoint.com) for complete details. Measure shoulder width to determine size. BPSFM and BPSL are approximately 2" wider in center than BPSFNM and BPSFNL.

BPSFM Anterior Trunk Support with Zip (13" - 16") .....	\$185.00 E0960	BPSFNM Anterior Trunk Support Narrow with Zip (13" - 16").....	\$185.00 E0960
BPSFL Anterior Trunk Support with Zip (16" - 19").....	\$185.00 E0960	BPSFNL Anterior Trunk Support Narrow with Zip (16" - 19").....	\$185.00 E0960

## MISCELLANEOUS ACCESSORIES

### BACK ACCESSORIES

ADJRGD Adjustable Width Back Rigidizing Bar.....	\$85.00 K0108	AASH Adj.Angle Stroller Handles (Outside Mount) .....	No Charge
1483 Non Adjustable Stroller Handle .....	\$90.00	AASHIN Adj.Angle Stroller Handles (Inside Mount) .....	No Charge

### MISCELLANEOUS ACCESSORIES

1209 Tire Pump.....	\$35.00 A9270	2218 IV Pole .....	\$175.00 K0105
P750 Oxygen Holder .....	\$150.00 E2208	Not available with LOSHEAR or LOSHEARWO.	
Not available with HDPG, LOSHEAR, LOSHEARWO, BH20, BH20B, ABH1720.		A1580 Rear Anti-Tippers.....	\$150.00
		BP02 Multi Pocket Backpack.....	\$85.00 A9270

## FRONT RIGGING OPTIONS

### FRONT RIGGING MOUNTING

STDMNT STD Front Rigging Mount .....	Standard
OFFMNT Offset Front Rigging Mount .....	No Charge
Not available with 12" - 16" seat depth range.	







FRONT RIGGINGS (PAIR)

Table listing various front rigging options such as 70LIFTNT 70° Lift Off Non-Taper Footrests, 80NTAPER 80° Swingaway Footrests, etc., with their respective prices and notes.

FOOTPLATE TYPE (PAIR)

Table listing footplate options such as 1651 Flip-Up Composite Footplates, A1562 Aluminum Footplates, A1564 Extra Large Aluminum Footplates, and AT5543 Adjustable Angle Flip-Up Footplates.

ADJUSTABLE ANGLE FOOTPLATE RIGHT (EACH) (BI-LATERAL CONTRACTURE FOOTREST ONLY)

Must select with BILAT.

Table listing adjustable angle footplate options for the right side, such as AAF325R, AAF425R, AAF525R, AAF625R, and AAF725R.

ADJUSTABLE ANGLE FOOTPLATE LEFT (EACH) (BI-LATERAL CONTRACTURE FOOTREST ONLY)

Must select with BILAT.

Table listing adjustable angle footplate options for the left side, such as AAF325L, AAF425L, AAF525L, AAF625L, and AAF725L.



Yes, you can.®

Knee-to-Heel Range (Measured from edge of seat pan to back of footplate)		
Front Rigging Type	Minimum	Maximum
Range adjusts in 1/2" increments.		
60°	13 1/2"	17 1/2"
60° MFX	10"	13 1/2"
*70°	13 1/2"	18 1/2"
*70° Tapered	13 1/4"	17.00"
80°	10.5"	16"
Bilateral Contracture STF125-STF16	5 1/2"	9"
Bilateral Contracture STF165-STF19	5 1/2"	13"
Range adjusts in 1" increments.		
70° MFX	6 1/2"	12 1/2"
80° MFX	4.5"	9"
90°	4 1/2"	7 1/2"
Infinite Range Adjustments		
*Elevating	14.00"	20.00"
Articulating	16.00"	24.00"
Contracture Platform STF125-STF16	7 1/4"	8 1/4"
Contracture Platform STF165-STF19	11 1/2"	16"

\* ALPT shifts range by 2"

**FOOTREST/LEGREST ACCESSORIES (PAIR)**

- 1337 Calf Strap ..... \$40.00 K0038
- 9342 "H" Strap ..... \$65.00 K0039
- FRMPRB ChairArmor Blue Neoprene Frame Protectors ..... \$90.00 A9270  
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- FRMPRR ChairArmor Red Neoprene Frame Protectors ..... \$90.00 A9270  
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- FRMPRBK ChairArmor Black Neoprene Frame Protectors ..... \$90.00 A9270  
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER

- ALPT Longer Pivot And Slide Tubes ..... \$84.00  
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- ALPT3 3" Extension For 90SWING Footrests ..... \$40.00 K0108
- SCLAMP Clamp-On Sector Blocks ..... \$75.00 K0108  
Not available with 70LIFTNT, 70LIFTMFX, 90LIFTMFX, AT5044LIFT, BILAT, 110CP. Adjustable angle. See Front Rigging Mounts image.
- CLFPDLRG Large Calf Pad for Elevating Leg Rests ..... No Charge  
Only available on Elevating leg rests.

**HEEL AND TOE LOOP (PAIR)**

- HLOOP Heel Loops ..... \$45.00 E0951
- 1345 Toe Loops ..... \$45.00 E0952



## SEATING

### SEATING CUSHION SELECTION

- ITSG Matrx Stabilite Gentle Contour Cushion ..... \$320.00  
E2605 E2606
- ITFG Matrx Flovair Gentle Contour Cushion ..... \$390.00  
E2607 E2608
- ITFM Matrx Flovair Max Contour Cushion ..... \$390.00 E2607  
E2608
- MA-VI Matrx Vi Cushion.....\$425.00 E2607 E2608

- PS-V Matrx PS Cushion..... \$345.00 E2605 E2606
- LC Matrx Libra Cushion 16"w - 20"w ..... \$485.00 E2624
- LC Matrx Libra Cushion 21"w ..... \$525.00 E2624
- LC Matrx Libra Cushion 22"w ..... \$595.00 E2625

### SEATING BACK SELECTION

- MPB-V Matrx PB Back .....\$625.00 E2615
- MPBD Matrx PB Deep Back ..... \$645.00 E2615
- PBE Matrx Elite Back..... \$625.00 E2620
- EDB Matrx Elite Deep Back..... \$760.00 E2620

- PBE-TR Matrx Elite TR..... \$845.00 E2620
- HDB Matrx PB Heavy Duty Back.....\$835.00 E2615 E2616
- DHDB Matrx PB Deep Heavy Duty Back ..... \$935.00 E2615 E2616

### HEADREST OPTIONS

- MEHSP-10 Matrx Elan Standard Pad 10" .....\$235.00 E0955

- MEHSP-14 Matrx Elan Standard Pad 14" .....\$235.00 E0955

### HEADREST MOUNTING HARDWARE

- MEMAH Elan Multi-Axis Removable Hardware .....\$220.00  
E1028

### FROGLEGS SUSPENSION FORKS

Front & Rear-Seat-to-Floor Heights Measured from floor to top of seat rail					
Rear Wheel	4" 4SPLT	5" 5IDC 5SPLT 5ALUM2	6" 6IDC 6ALUM2	6" 1259	8" 8IDC 8I5DC
12"	13.50"	15.50"	15.50" 16.00"	16.50" 17.00" 17.50" 18.00"	N/A
16"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00"	N/A	17.00" 17.50" 18.00"	18.50" 19.00"
18"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00" 15.50"	15.50" 16.00"	18.00"	18.50" 19.00"
20"	15.00"	15.00" 15.50"	15.50" 16.00"	16.50" 17.00" 17.50"	N/A
22"	N/A	N/A	16.00"	16.50" 17.00" 17.50" 18.00"	18.50"
24"	N/A	N/A	N/A	16.50" 17.00" 17.50" 18.00"	18.50" 19.00"

#### Basic Chair Includes:

- Solara Tilt-in-Space Frame
- Seat-to-Floor Height Selection
- Tilt Option (FRLT, FRRT, THR)
- Casters (4SPLT, 5SPLT, 5IDC, 6IDC, 8IDC)
- Frame Width (12"-19")
- Wheels
- Axles (A1577, PAXLE)
- Urethane Treaded Tires
- Handrims (AHDRM, CHDRM)
- Center-of-Gravity Setting - Based on seat depth and seat depth range. See page 7 for details.
- Frame Finish
- Two-Point Conventional Fixed Arms (FXDESK or FXFULL)
- Rear Anti-Tippers
- Aluminum or Steel Reinforced Back Canes
- Wheel Locks (PULLLOCK, PUSHLOCK)
- Back Height Selection
- Front Riggings (SWING1319, SWING914, 70NTAPER, 70MFX, 70STAPER, 90SWING)
- Fold Down Back
- Flip-up Composite Footplates

Specifications	
Frame Style	Tilt-in-Space
Overall Width	Seat width + 8 1/4" (12" wheels) Seat width + 9 1/4" (18-24" wheels)
Seat Widths	12" - 24"
Seat Depths	12" - 22"
Back Height	17" - 20", 20", 24"
Adjustable Angle Back	0°, 5°, 10°, 15°, 20°, 25°, 30°
Tilt Range	-5° to 50°
Recliner Back	160°
Weight Limit	300 lbs./400 lbs. with Heavy-Duty Package
Weight	34 1/2 lbs*
Warranty	Lifetime on frame and crossbraces

\*Base chair without Accessories

#### Available Through Our Specials Departments

- 20" High Reclining Back SM-050503023551
- 80° Footrest Hanger SM-200411121756
- Front Rigging Mount: C-090731110749  
Standard on one side and offset on the opposite side

#### Tilt Assist Explained

The Solara® 3G wheelchair is shipped from the factory pre-configured with one of three "Tilt Assist" settings: Front, Mid or Rear, based solely on the seat depth range and seat depth ordered. This setting describes how the user's weight is distributed on the frame itself, i.e., does one front load the chair, rear load or is the weight evenly distributed. As compared to the 2G, the 3G requires up to 30% less force needed to tilt because of these settings. If factors other than seat depth or seat depth range come into play that affect the weight distribution of the chair, the provider has the ability to change the Tilt Assist setting in the field by simply adjusting 2 screws on each side of the frame.

Tilt Assist Factory Setting		
	Seat Depth Range 15" - 19"	Seat Depth Range 18" - 22"
D15	Front Setting	
D16	Front Setting	
D17	Mid Setting	
D18	Rear Setting	Front Setting
D19	Rear Setting	Front Setting
D20		Mid Setting
D21		Rear Setting
D22		Rear Setting


The factory setting may need to be adjusted if a seating system is added. The factory setting does not account for conditions outside of the seat depth and seat depth range which may cause more weight to be distributed to the front or back of the chair.

When adjusting the Tilt Assist there are 4 bolts to change (2 on each side.) See the illustration to the left and chart above.

The Rear setting has the rocker tube telescoped entirely into the seat frame support tube where the adjustment bolt is in the bottom hole. Along the base frame the bolt is in the top hole. (e.g. accommodates a seat depth that is in the highest end of the seat depth range)

For the Mid setting the adjustment bolt is simply in the middle hole along the rocker tube which telescopes out and the middle hole along the base frame. (e.g. optimal for a seat depth in the middle of the seat depth range)

The Front setting has the bolt in the top hole of the rocker tube which is fully telescoped out of the seat frame support tube and in the bottom hole along the base frame. (e.g. accommodates a seat depth in the lowest end of the seat depth range)



## Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.