




Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

Thank you for using Invacare eForms.

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit www.invacare.com/eForms



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-678-4682



Account Information

Request Type: Quote Order

Date: _____

Account #: _____

Company: _____

SHIP TO
Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Country: _____

Purchase Order #: _____

CONTACT
Name: _____

Back Up Contact: _____

Phone: _____

Email: _____

Comments: _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTIONS

FRAME TYPE

SOLARA3G Solara 3G Tilt-in-Space wheelchair\$ 3,228.00 E1161

A seating system must be purchased. No seating system is provided. 300 lb standard weight capacity. See www.invacare.com/Solara3G for additional information when configuring this chair.

CUSTOM PACKAGES

HDPG Heavy Duty Package (400 lbs. Weight Capacity) \$420.00 K0108

Reinforced axles and crossbraces. Weight limit of 400 lbs.

SUSPENSION SELECTION

SUSPENSION

DROCK Seating Dynamics Rocker Back.....\$2,190.00 E2398

Not available with HDPG, LOSHEAR, LOSHEARWO, RCANE, FLDBACK. See www.seatingdynamics.com for more details.

FLSF Frog Legs Suspension Forks.....\$439.00 E1015

Not available with HDPG. See chart page 7.

SPACK Suspension Package\$1,599.00 E1015

Not available with HDPG, LOSHEAR, LOSHEARWO, RCANE, FLDBACK. Includes DROCK and FLSF.

Transport Tie Down

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tiedown systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

TRANSPORT/VENT TRAY SELECTION

TRANSPORT TIE DOWN

TRBKTS Transport Brackets \$331.00 K0108

Not available with AVENT, BACKVENT, P750, 2218. Now available as a service part kit for existing models.

VENT TRAY OPTION

AVENT Ventilator Tray (each).....\$326.00 E1029

Not available with HDPG. Available with 16" wide frame and above and 18" Seat-to-Floor Height and above.

ATBOX Ventilator Battery Box (each)\$277.00

BACKVENT Laptop Vent Tray \$1,698.00 E1030

Not available with HDPG, LOSHEAR, P750. Available with 15" wide frame and above.



FRAME SELECTIONS

TILTING OPTIONS

- FRLT Cable-Free Foot Release-Left Mount..... No Charge
- FRRT Cable-Free Foot Release-Right Mount.....No Charge
- THR Traditional Hand Release Standard

- TLL Tilt Lock Levers..... No Charge
- TLC Tilt Lock Stops (pair)..... \$39.00
Additional pair of tilt lock stops, one pair comes standard with chair.

SEAT WIDTH RANGE

- SWR1216 12" – 16" Seat Width Range No Charge
Not available with HDPG.
- SWR1418 14" – 18" Seat Width Range..... No Charge
- SWR1620 16" – 20" Seat Width Range No Charge

- SWR1822 18" – 22" Seat Width Range No Charge
- SWR2024 20" – 24" Seat Width Range..... No Charge

SEAT WIDTH

- W12 12" Wide No Charge
- W13 13" Wide No Charge
- W14 14" Wide No Charge
- W15 15" Wide No Charge
- W16 16" Wide No Charge
- W17 17" Wide No Charge
- W18 18" Wide No Charge

- W19 19" Wide No Charge
- W20 20" Wide \$552.00 E2201
- W21 21" Wide\$552.00 E2201 K0108
- W22 22" Wide\$552.00 E2201 K0108
- W23 23" Wide\$552.00 E2201 K0108
- W24 24" Wide\$552.00 E2202 K0108

SEAT DEPTH RANGE

- SDR1216 12" – 16" Seat Depth Range..... No Charge
Tilt Assist feature not available on this frame.
- SDR1519 15" – 19" Seat Depth Range No Charge

- SDR1822 18" – 22" Seat Depth Range..... No Charge

SEAT DEPTH

- D12 12" Seat Depth No Charge
- D13 13" Seat Depth No Charge
- D14 14" Seat Depth No Charge
- D15 15" Seat Depth No Charge
- D16 16" Seat Depth No Charge
- D17 17" Seat Depth No Charge

- D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge
- D20 20" Seat Depth..... \$607.00 E2203
- D21 21" Seat Depth.....\$607.00 E2203 K0108
- D22 22" Seat Depth\$1,092.00 E2204 K0108



Yes, you can.

FRAME FINISH

FRAME COLORS

- 24P Wet Black No Charge
- 29P White Snow No Charge
- 60P Silver Metallic No Charge
- 115P Black Prism No Charge
- 105P Midnight Blue No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 30P Sunny Yellow No Charge

- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge
- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge
- Touch Up Paint \$23.00

Touch up paint is not available in all colors.

Touch Up Color

Discontinued frame colors are available through our specials department. C-070109030544

CASTER, WHEEL AND TIRE SELECTIONS

SEAT-TO-FLOOR HEIGHT

- STF125 Seat-To-Floor 12 1/2" (10 1/2" w/DPSEAT)..... No Charge
- STF13 Seat-To-Floor 13" (11" w/DPSEAT)..... No Charge
- STF135 Seat-To-Floor 13 1/2" (11 1/2" w/DPSEAT) No Charge
- STF14 Seat-To-Floor 14" (12" w/DPSEAT) No Charge
- STF145 Seat-To-Floor 14 1/2" (12 1/2" w/DPSEAT)..... No Charge
- STF15 Seat-To-Floor 15" (13" w/DPSEAT) No Charge
- STF155 Seat-To-Floor 15 1/2" (13 1/2" w/DPSEAT)..... No Charge

- STF16 Seat-To-Floor 16" (14" w/DPSEAT) No Charge
- STF165 Seat-To-Floor 16 1/2" (14 1/2" w/DPSEAT)..... No Charge
- STF17 Seat-To-Floor 17" (15" w/DPSEAT) No Charge
- STF175 Seat-To-Floor 17 1/2" (15 1/2" w/DPSEAT)..... No Charge
- STF18 Seat-To-Floor 18" (16" w/DPSEAT) No Charge
- STF185 Seat-To-Floor 18 1/2" (16 1/2" w/DPSEAT) No Charge
- STF19 Seat-To-Floor 19" (17" w/DPSEAT) No Charge

CASTER OPTIONS (PAIR)

- 4SPLT 4" x 1" Composite Urethane Light Up..... No Charge
Not available with HDPG.
- 4ALUM2 4" x 1.4" Aluminum Soft Roll \$139.00
- 5SPLT 5" x 1" Composite Urethane Light Up..... No Charge
Not available with HDPG.
- 51DC 5" x 1" Urethane Casters..... No Charge
Not available with HDPG.
- 5ALUM2 5" x 1.4" Aluminum Soft Roll \$139.00

- 61DC 6" x 1" Urethane Caster..... No Charge
Not available with HDPG.
- 1259 6" x 2" Semi-Pneumatic Casters..... \$95.00
Not available with HDPG.
- 6ALUM2 6" x 1.4" Aluminum Soft Roll \$139.00
- 81DC 8" x 1" Urethane Casters No Charge
- 815DC 8" x 1.5" Urethane Caster No Charge



Yes, you can.®

WHEEL OPTIONS (PAIR)

- 12COM 12" Composite Wheels No Charge
- 16COM 16" Composite Wheels No Charge
Not available with urethane tires. Only available with aluminum handrims.
- 18COM 18" Composite Wheels No Charge
Not available with urethane tires. Only available with aluminum handrims.

- 20COM 20" Composite Wheels No Charge
- 22COM 22" Composite Wheels No Charge
- 24COM 24" Composite Wheels No Charge

AXLE TYPE (PAIR)

A1577 Quick Release Axles No Charge

PAXLE Permanent Axles Standard

TIRE OPTIONS (PAIR)

- PNEU Pneumatic Tires \$95.00
Not available with HDPG.
- PNEUFF Flat Free Inserts \$95.00 E2213
- URETH Urethane (Treaded) Tires No Charge

HANDRIM OPTIONS (PAIR)

- CHDRM Chrome Plated Handrims No Charge
- 1512 Plastic Coated Aluminum Handrims \$128.00
- AT-1225-8 Plastic Coated with 8 Vertical Projections \$326.00
.....
Only available with 24COM7.
- AT-1224-8 Plastic Coated with 8 Oblique Projections \$326.00
.....
Only available with 24COM7.
- AT1225-12 Plastic Coated with 12 Vertical Projections \$326.00
.....
Only available with 24COM7.
- L8702 Omit Handrims \$0.00

Front and Rear Seat-to-Floor Heights (measured from Floor to Top of Seat Rail)					
Rear Wheel	4" 4SPLT 4ALUM2	5" 51DC 5ALUM2 5SPLT	6" 61DC 6ALUM2	6" 1259	8" 81DC 815DC
12" Pneu PneuFF	12 1/2" 13" 13 1/2" 15 1/2" 16"	13 1/2" 15 1/2" 16" 16 1/2"	15 1/2" 16" 16 1/2" 17"	15 1/2" 16" 16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18"
12" Urethane	12 1/2" 14 1/2" 15" 15 1/2" 16" 16"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	14 1/2" 15" 15 1/2" 16" 16 1/2"	16 1/2" 17"
16"	12 1/2" 13" 13 1/2" 14" 14 1/2" 15"	13 1/2" 14" 14 1/2" 15" 15"	14 1/2" 15" 17"	14 1/2" 15" 17" 17 1/2" 18"	17" 17 1/2" 18" 18 1/2" 19"
18"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	13 1/2" 14" 14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16"	14 1/2" 15" 15 1/2" 16" 18"	18" 18 1/2" 19"
20"	15" 15 1/2" 16"	15" 15 1/2" 16" 16 1/2"	15" 15 1/2" 16" 16 1/2" 17"	15" 15 1/2" 16" 16 1/2" 17" 17 1/2"	16 1/2" 17" 17 1/2"
22"	16"	16" 16 1/2"	16" 16 1/2" 17"	16" 16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18" 18 1/2"
24"	N/A	16 1/2"	16 1/2" 17"	16 1/2" 17" 17 1/2" 18"	16 1/2" 17" 17 1/2" 18" 18 1/2" 19"

All Dimensions Are Approximate



Yes, you can.®

WHEEL LOCKS

FOOTBRK Foot Operated Push To Lock Wheel Lock\$83.00 K0108
 Not available with AVENT.

PULLLOCK Pull To Lock Wheel Locks (pair) No Charge

PUSHLOCK Push To Lock Wheel Locks (pair).....No Charge

EXTLOCK Push To Lock with 6 Bolt-On Ext (pair).....\$83.00
 E0961

FOOTHUB Foot Operated Hublock \$221.00 K0108
 Not available with 16COM.

ATHUB Cane Mount Hub Lock \$221.00 K0108
 Not available with LOSHEAR or LOSHEARWO, 16COM.

USEHUB Seat Rail Mount Hub Lock \$221.00 K0108
 Only available with TAD of TAF. Not available with 16COM.

ARM AND BACK SELECTIONS

SEAT STYLE

FLSEAT or DPSEAT are solid seat pans and do not come standard with chair. Both pans grow in width and require no new hardware or parts when converting from flush to drop or vice versa.

FLSEAT Flush Seat Pan\$177.00 E2231

DPSEAT Drop Seat Pan \$177.00 E2231
 Accommodates 2" lower seat to floor heights. Use a cushion not less than 3" thick.

NOSEAT No Seat PanNo Charge

Seat Pan Growth Range (1" increments)	
Chair Width	Seat Pan Width Range
W12	12"-14"
W13	12"-14"
W14	13"-16"
W15	13"-16"
W16	15"-18"
W17	15"-18"
W18	17"-20"
W19	17"-20"
W20	19"-22"
W21	19"-22"
W22	21"-24"
W23	21"-24"
W24	23"-24"

Arm Heights (measured from Top of Seat Rail to Top of Armpad)		
Arm Type	Low	High
TAD	6 1/2" - 12 1/2"	9" - 17"
TAF	1/2" increments	1" increments
FXDESK	9 1/2"	12"
FXFULL		
ADJDESK	9 1/2" - 13 1/2"	12" - 16"
ADJFULL	1" increments	1" increments
FADJDESK	9" - 13"	12" - 16"
FADJFULL	1" increments	1" increments

ARM HEIGHT

SARMHT Low Arm Height..... No Charge

HARMHT High Arm Height..... No Charge

CANTILEVER ARMS

Not available with LOSHEAR or LOSHEARWO

AT99N Cantilever Non Locking Bolt on ArmsNo Charge

AT55D Cantilever Locking Desk Length Bolt on Arms \$293.00

AT55F Cantilever Locking Full Length Bolt on Arms..... \$293.00

T-ARMS (PAIR)

TAD Adjustable Height, Desk Length\$315.00 E0973

TAF Adjustable Height, Full Length\$315.00 E0973

TWO-POINT, FIXED HEIGHT ARMS (PAIR)

Available on 15" seat depth and above.

FXDESK Conventional, Desk LengthNo Charge

FXFULL Conventional, Full LengthNo Charge



Yes, you can.

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

TWO-POINT, ADJUSTABLE HEIGHT ARMS (PAIR)

Available on 15" seat depth and above.

ADJDESK Conventional, Desk Length \$293.00 E0973

ADJFULL Conventional, Full Length \$293.00 E0973

FADJDESK Flip-Back, Desk Length \$326.00 E0973

Not available with LOSHEAR, LOSHEARWO or HDPG.

FADJFULL Flip-Back, Full Length \$326.00 E0973

Not available with LOSHEAR, LOSHEARWO or HDPG.

BACK TYPE

SCANE Standard Aluminum Back Cane Standard

RCANE Steel-Reinforced Back Cane No Charge

LOSHEAR Low Shear Reclining Back with Upholstery
..... \$1,004.00 E1226

Not available with HDPG. Comes with back upholstery and AASH. Available on 14"-20" seat width.

Not available with Foot Operated Wheel Lock

LOSHEARWO Low Shear Reclining Back without Upholstery
..... \$971.00 E1226

Not available with HDPG, headrest and pillow. Comes standard with AASH.

Not available with Foot Operated Wheel Lock

RECEXT Reclining Headrest Pillow and Extension
..... \$95.00 E0966

RLL Recliner Lockout Levers No Charge

Not available with SCANE or RCANE.

RLC Recliner Lock Stops \$39.00

Back Cane Heights (measured from Top of Seat Rail)		
Back Cane Type	Inches to 10 degree bend	Total Height
BH20	Straight Cane	20.5"
BH20B/BH24B	15 1/2"	24.5"
BH24	Straight Cane	24 1/2"
LOSHEAR	Straight Cane	26"
LOSHEAR + Headrest	Straight Cane	36"
ABH1720	16.2" - 21.2"	29.3" - 34.3"

BACK HEIGHT

BH20B Back Height 20" with bend No Charge
Not available with P750.

BH20 Strt Cane w Push Handle Back Ht 20" No Charge
A9900
Not available with P750.

BH24B Back Height 24" with bend No Charge

BH24 Back Height 24" No Charge

ABH1720 Adj Height Back Canes (17" – 20") No Charge
Not available with RCANE (because it consists of heavy duty construction), 1483, ATHUB, P750.

BACK CANE OPTIONS

FXBACK Non-Fold Down Back Canes Standard

FLDBACK Fold-Down Back Canes No Charge
Not available with LOSHEAR, DROCK.

POSITIONING STRAPS (EACH)

1515 Push Button Buckle Seat Positioning Strap \$63.00
E0978

1321BK Hook and Loop Seat Positioning Strap \$63.00
E0978

9515 Push Button Chest Positioning Strap \$123.00 E0960
Not available with LOSHEAR or LOSHEARWO.

7321BK Hook and Loop Chest Positioning Strap \$130.00
E0960
Not available with LOSHEAR or LOSHEARWO.



Yes, you can.®

BODYPOINT® PADDED CALF STRAPS

See www.bodypoint.com for complete details. Measure chair width to determine size.

BPCALFS Padded Calf Straps (13" - 15")..... \$76.00 K0038

BPCALFL Padded Calf Straps (18" - 22") \$76.00 K0038

BPCALFM Padded Calf Straps (15" - 18") \$76.00 K0038

BODYPOINT® ANKLE HUGGERS®

See www.bodypoint.com for complete details. Measure circumference above ankle to determine size.

BPAHGS Side - Release Ankle Huggers (6.5" - 8").....\$199.00
K0108

BPAHGL Side - Release Ankle Huggers (8.5" - 11.5").....\$199.00
K0108

BPAHGM Side - Release Ankle Huggers (7.5" - 9").....\$199.00
K0108

BPAHGXL Side - Release Ankle Huggers (11" - 13")\$199.00
K0108

MISCELLANEOUS ACCESSORIES

BACK ACCESSORIES

ADJRGD Adjustable Width Back Rigidizing Bar.....\$95.00 K0108

AASH Adj.Angle Stroller Handles (Outside Mount)
.....No Charge

1483 Non Adjustable Stroller Handle\$100.00

AASHIN Adj.Angle Stroller Handles (Inside Mount)
.....No Charge

MISCELLANEOUS ACCESSORIES

1209 Tire Pump.....\$39.00 A9270

2218 IV Pole.....\$194.00 K0105

Not available with LOSHEAR or LOSHEARWO.

P750 Oxygen Holder\$166.00 E2208

A1580 Rear Anti-Tippers.....\$166.00

Not available with HDPG, LOSHEAR, LOSHEARWO, BH20, BH20B, ABH1720.

BP02 Multi Pocket Backpack.....\$95.00 A9270

FRONT RIGGING OPTIONS

FRONT RIGGING MOUNTING

STDMNT STD Front Rigging MountStandard

OFFMNT Offset Front Rigging Mount No Charge

Not available with 12" - 16" seat depth range.





Yes, you can.®

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

FRONT RIGGINGS (PAIR)

70LIFTNT 70° Lift Off Non-Taper Footrests.....	\$105.00
70LIFTMFX 70° Lift Off MFX Footrests.....	\$105.00
Must select AT5543.	
90LIFTMFX 90° Lift Off MFX Footrests	\$105.00
Must select AT5543.	
SWING1319 60° Swingaway Footrests.....	No Charge
SWING914 60° Swingaway Footrests MFX.....	No Charge
E0990	
70NTAPER 70° Swingaway Footrests	No Charge
70MFX 70° Swingaway Footrests MFX	No Charge E0990
Must select AT5543.	
70STAPER 70° Tapered Footrests	No Charge E0990
Available on seat widths 16" and above.	

80NTAPER 80° Swingaway Footrests.....	No Charge
80MFX 80° Swingaway Footrests MFX.....	No Charge
Must select AT5543.	
90SWING 90° Swingaway Footrests.....	No Charge
Must select AT5543.	
AT5044LIFT Lift Off Elevating Legrests	\$338.00 E0990
AT5044 Elevating Legrests	\$338.00 E0990
BILAT Bi-Lateral Contracture Footrests.....	\$387.00
Not available with HDPG OR 12 1/2" STF Height.	
ATL4A Center Pivot Articulating Legrests	\$387.00 K0053
110CP Contracture Platform	\$387.00
Not available with HDPG. Available on 16"-22" seat widths.	

FOOTPLATE TYPE (PAIR)

1651 Flip-Up Composite Footplates	No Charge
Not available with Bodypoint Ankle Huggers.	
A1562 Aluminum Footplates	\$94.00
Not available with Bodypoint Ankle Huggers.	

A1564 Extra Large Aluminum Footplates	\$94.00 K0041
Not available with Bodypoint Ankle Huggers.	
AT5543 Adjustable Angle Flip-Up Footplates	\$200.00
K0040	

ADJUSTABLE ANGLE FOOTPLATE RIGHT (EACH) (BI-LATERAL CONTRACTURE FOOTREST ONLY)

Must select with BILAT.

K0040

AAF325R Adj. Angle Footplate Right 3 1/4" W	\$105.00
K0040	
AAF425R Adj. Angle Footplate Right 4 1/4" W	\$105.00
K0040	
AAF525R Adj. Angle Footplate Right 5 1/4" W	\$105.00

AAF625R Adj. Angle Footplate Right 6 1/4" W	\$105.00
K0040	
AAF725R Adj. Angle Footplate Right 7 1/4" W	\$105.00
K0040	

ADJUSTABLE ANGLE FOOTPLATE LEFT (EACH) (BI-LATERAL CONTRACTURE FOOTREST ONLY)

Must select with BILAT.

AAF325L Adj. Angle Footplate Left 3 1/4" W	\$105.00 K0040
AAF425L Adj. Angle Footplate Left 4 1/4" W	\$105.00 K0040
AAF525L Adj. Angle Footplate Left 5 1/4" W	\$105.00 K0040

AAF625L Adj. Angle Footplate Left 6 1/4" W	\$105.00 K0040
AAF725L Adj. Angle Footplate Left 7 1/4" W	\$105.00 K0040



Yes, you can.®

Knee-to-Heel Range (Measured from edge of seat pan to back of footplate)		
Front Rigging Type	Minimum	Maximum
Range adjusts in 1/2" increments.		
60°	13 1/2"	17 1/2"
60° MFX	10"	13 1/2"
*70°	13 1/2"	18 1/2"
*70° Tapered	13 1/4"	17.00"
80°	10.5"	16"
Bilateral Contracture STF125-STF16	5 1/2"	9"
Bilateral Contracture STF165-STF19	5 1/2"	13"
Range adjusts in 1" increments.		
70° MFX	6 1/2"	12 1/2"
80° MFX	4.5"	9"
90°	4 1/2"	7 1/2"
Infinite Range Adjustments		
*Elevating	14.00"	20.00"
Articulating	16.00"	24.00"
Contracture Platform STF125-STF16	7 1/4"	8 1/4"
Contracture Platform STF165-STF19	11 1/2"	16"

* ALPT shifts range by 2"

FOOTREST/LEGREST ACCESSORIES (PAIR)

- 1337 Calf Strap..... \$45.00 K0038
- 9342 "H" Strap..... \$73.00 K0039
- FRMPRB ChairArmor Blue Neoprene Frame Protectors \$100.00 A9270
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- FRMPRR ChairArmor Red Neoprene Frame Protectors \$100.00 A9270
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER
- FRMPRBK ChairArmor Black Neoprene Frame Protectors \$100.00 A9270
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER

ALPT Longer Pivot And Slide Tubes..... \$89.00
Only available with 70LIFTNT, SWING1319, 70NTAPER, 70STAPER

ALPT3 3" Extension For 90SWING Footrests..... \$45.00
K0108

SCLAMP Clamp-On Sector Blocks\$83.00 K0108

Not available with 70LIFTNT, 70LIFTMFX, 90LIFTMFX, AT5044LIFT, BILAT, 110CP. Adjustable angle. See Front Rigging Mounts image.

CLFPDLRG Large Calf Pad for Elevating Leg Rests..... No Charge

Only available on Elevating leg rests.

HEEL AND TOE LOOP (PAIR)

HLOOP Heel Loops\$51.00 E0951

1345 Toe Loops\$51.00 E0952



Yes, you can.®

SEATING

SEATING CUSHION SELECTION

MA-VI Matrix Vi Cushion.....\$470.00 E2607 E2608
 PS-V Matrix PS Cushion.....\$382.00 E2605 E2606
 LC Matrix Libra Cushion 16"w - 20"w..... \$536.00 E2624

LC Matrix Libra Cushion 21"w or 22"w.....\$536 E2624

Must select appropriate 21" or 22" up-charge

LC Matrix Libra Cushion 21"w up-charge..... \$45.00 E2624

LC Matrix Libra Cushion 22"w up-charge..... \$116.00 E2625

SEATING BACK SELECTION

MPB-V Matrix PB Back\$690.00 E2615
 MPBD Matrix PB Deep Back \$712.00 E2615
 PBE Matrix Elite Back..... \$690.00 E2620

EDB Matrix Elite Deep Back..... \$838.00 E2620

HDB Matrix PB Heavy Duty Back.....\$921.00 E2615 E2616

DHDB Matrix PB Deep Heavy Duty Back \$1,032.00 E2615 E2616

HEADREST OPTIONS

MEHSP-10 Matrix Elan Standard Pad 10"\$260.00 E0955

MEHSP-14 Matrix Elan Standard Pad 14"\$260.00 E0955

HEADREST MOUNTING HARDWARE

MEMAH Elan Multi-Axis Removable Hardware\$243.00
 E1028

FROGLEGS SUSPENSION FORKS

Front & Rear-Seat-to-Floor Heights Measured from floor to top of seat rail					
Rear Wheel	4" 4SPLT	5" 5IDC 5SPLT 5ALUM2	6" 6IDC 6ALUM2	6" 1259	8" 8IDC 8I5DC
12"	13.50"	15.50"	15.50" 16.00"	16.50" 17.00" 17.50" 18.00"	N/A
16"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00"	N/A	17.00" 17.50" 18.00"	18.50" 19.00"
18"	13.50" 14.00" 14.50" 15.00"	14.50" 15.00" 15.50"	15.50" 16.00"	18.00"	18.50" 19.00"
20"	15.00"	15.00" 15.50"	15.50" 16.00"	16.50" 17.00" 17.50"	N/A
22"	N/A	N/A	16.00"	16.50" 17.00" 17.50" 18.00"	18.50"
24"	N/A	N/A	N/A	16.50" 17.00" 17.50" 18.00"	18.50" 19.00"

Basic Chair Includes:

- Solara Tilt-in-Space Frame
- Seat-to-Floor Height Selection
- Tilt Option (FRLT, FRRT, THR)
- Casters (4SPLT, 5SPLT, 5IDC, 6IDC, 8IDC)
- Frame Width (12"-19")
- Wheels
- Axles (A1577, PAXLE)
- Urethane Treaded Tires
- Handrims (AHDRM, CHDRM)
- Center-of-Gravity Setting - Based on seat depth and seat depth range. See page 7 for details.
- Frame Finish
- Two-Point Conventional Fixed Arms (FXDESK or FXFULL)
- Rear Anti-Tippers
- Aluminum or Steel Reinforced Back Canes
- Wheel Locks (PULLLOCK, PUSHLOCK)
- Back Height Selection
- Front Riggings (SWING1319, SWING914, 70NTAPER, 70MFX, 70STAPER, 90SWING)
- Fold Down Back
- Flip-up Composite Footplates



Yes, you can.®

Invacare® Solara™ 3G - Adult Tilt-in-Space Wheelchair

Specifications	
Frame Style	Tilt-in-Space
Overall Width	Seat width + 8 1/4" (12" wheels) Seat width + 9 1/4" (18-24" wheels)
Seat Widths	12" - 24"
Seat Depths	12" - 22"
Back Height	17" - 20", 20", 24"
Adjustable Angle Back	0°, 5°, 10°, 15°, 20°, 25°, 30°
Tilt Range	-5° to 50°
Recliner Back	160°
Weight Limit	300 lbs./400 lbs. with Heavy-Duty Package
Weight	34 1/2 lbs*
Warranty	Lifetime on frame and crossbraces

*Base chair without Accessories

Available Through Our Specials Departments

- 20" High Reclining Back SM-050503023551
- 80° Footrest Hanger SM-200411121756
- Front Rigging Mount: C-090731110749
Standard on one side and offset on the opposite side

Tilt Assist Explained

The Solara® 3G wheelchair is shipped from the factory pre-configured with one of three "Tilt Assist" settings: Front, Mid or Rear, based solely on the seat depth range and seat depth ordered. This setting describes how the user's weight is distributed on the frame itself, i.e., does one front load the chair, rear load or is the weight evenly distributed. As compared to the 2G, the 3G requires up to 30% less force needed to tilt because of these settings. If factors other than seat depth or seat depth range come into play that affect the weight distribution of the chair, the provider has the ability to change the Tilt Assist setting in the field by simply adjusting 2 screws on each side of the frame.

Tilt Assist Factory Setting		
	Seat Depth Range 15" - 19"	Seat Depth Range 18" - 22"
D15	Front Setting	
D16	Front Setting	
D17	Mid Setting	
D18	Rear Setting	Front Setting
D19	Rear Setting	Front Setting
D20		Mid Setting
D21		Rear Setting
D22		Rear Setting


The factory setting may need to be adjusted if a seating system is added. The factory setting does not account for conditions outside of the seat depth and seat depth range which may cause more weight to be distributed to the front or back of the chair.

When adjusting the Tilt Assist there are 4 bolts to change (2 on each side.) See the illustration to the left and chart above.

The Rear setting has the rocker tube telescoped entirely into the seat frame support tube where the adjustment bolt is in the bottom hole. Along the base frame the bolt is in the top hole. (e.g. accommodates a seat depth that is in the highest end of the seat depth range)

For the Mid setting the adjustment bolt is simply in the middle hole along the rocker tube which telescopes out and the middle hole along the base frame. (e.g. optimal for a seat depth in the middle of the seat depth range)

The Front setting has the bolt in the top hole of the rocker tube which is fully telescoped out of the seat frame support tube and in the bottom hole along the base frame. (e.g. accommodates a seat depth in the lowest end of the seat depth range)



Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.