



Yes, you can.®

### Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

#### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

#### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® TDX® SR Wheelchair with Captain's Seat

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### CHAIR TYPE

- TDXSRV TDX SR with MK6i™ ..... \$9,950.00 K0869  
**TDXSV comes standard with: 300 lb. user weight capacity, 7 MPH top speed, GB TrueTrack motor package, MK6TT Controller and Quiet Stability Lock. Base width 25 1/2"**
- TDXSRV-HD TDX SR Heavy Duty with MK6i ..... \$10,500.00  
**TDXSRV-HD comes standard with: 400 lb. user weight capacity, 7 MPH top speed, GB TrueTrack motor package, MK6TT Controller and Quiet Stability Lock. Base width 25 1/2"**
- TDXSRSEAT TDXSR Seat Only ..... \$755.00

### USER WEIGHT

U400 User Weight <400 lbs. .... \$400.00

### WHEEL OPTIONS (PAIR)

1431-3 14" x 3" with Foam Filled Inserts ..... Standard

### FORK OPTION

FKPKG One-Sided Fork Package ..... Standard

### FRAME FINISH

#### CONTEMPORARY

- 24P Wet Black ..... No Charge
- 60P Silver Metallic ..... No Charge
- 105P Midnight Blue ..... No Charge
- 115P Black Prism ..... No Charge

#### TOO HOT

- 30P Sunny Yellow ..... No Charge
- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge

#### WAY COOL

- 62P Electric Blue ..... No Charge
- 127P Grasshopper Green ..... No Charge
- 116P Sunkissed Orange ..... No Charge
- 117P Tangerine Red ..... No Charge
- 29P White Snow ..... No Charge

Discontinued colors 71P Silver Vein, 104P Electric Teal, 121P Lolly Pop Blue, 122P Cosmic Blue, and 125P Emerald Green can be order through specials number C-070109031327

Group 3 - K0869 Basic Package	
TDX SR with MK6i and Van Seat	
400 lbs. Weight Capacity	\$400.00
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
TrueTrack Gearless Brushless GB Motors	Standard
Flip Back Fixed Height Arm	Standard
MK690ACC TrueTrack Controller	Standard
SPJPLUS SPJ™+ Single Drive Joystick	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard

### TRANSPORT TIE DOWN

TRBKTS Wheelchair Transport Brackets ..... \$250.00

#### TRANSPORT TIE-DOWN\*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

**NOTE: All specifications and dimensions are approximate.  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® TDX® SR Wheelchair with Captain's Seat

## VAN SEAT (STANDARD DIMENSIONS: 18"W x 16"D)

- VW16 16" Wide Van Seat ..... No Charge
  - VD16 16" Depth ..... No Charge
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW18 18" Wide Van Seat ..... Standard
  - VD16 16" Depth ..... Standard
  - VD17 17" Depth ..... No Charge
  - VD18 18" Depth ..... No Charge
- VW20 20" Wide Van Seat ..... \$329.00
  - VD18 18" Seat Depth ..... No Charge
  - VD19 19" Seat Depth ..... No Charge
  - VD20 20" Seat Depth ..... \$529.00
- VW22 22" Wide Van Seat ..... \$329.00
  - VD18 18" Seat Depth ..... No Charge
  - VD19 19" Seat Depth ..... No Charge
  - VD20 20" Seat Depth ..... \$529.00

## VAN UPHOLSTERY

- GV Grey Vinyl (Van Seat) ..... No Charge

## SEAT STYLE

- Clear
- SEMIR Semi Recline ..... Standard
- SEMIRB Semi Recline with Solid Base ..... No Charge

## ARM TYPE

Arms are flip-up and adjustable.

- Clear
- VAH Full Length, Adjustable Height ..... \$230.00 E0973
- VFH Full Length, Fixed Height ..... Standard
- VAHD Desk Length, Adjustable Height ..... \$230.00 E0973
- VFHD Desk Length, Fixed Height ..... No Charge

## SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap ..... \$50.00
- 1515 Push Button Style Seat Positioning Strap ..... Standard

## MK6i® CONTROLLER

- MK6TT TrueTrack Controller ..... Standard

## EXPANDABLE CONTROLLERS SYSTEMS

- Clear
- PWH Harness req. for Expandable Electronics ..... \$350.00 E2313  
Required to interface joystick or MK6i display with expandable controller.
- MPJM6U MPJ+ Multiple Drive Joystick ..... \$750.00 E2377

## JOYSTICK MOUNTING HARDWARE

- QLM6 Quad Link Retract. Jystk Mnt CMPJ+ Van ..... \$230.00 E1028
- ARM250 Stealth Height Adj Joystick Mount ..... \$210.00 E1028
- LEFT Left Handed Mounted Joystick ..... No Charge

## ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card ..... \$475.00
- 1813SD Pro Memory Card (USB Ready) ..... \$90.00
- 1560 "T" Handle Flexible Joystick Ext. .... \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. .... \$77.00 E2323
- A24VPS 24V MK6i Auxiliary Power Source ..... \$200.00
- MPY Monoport "Y" Cable+ or CMPJ+ or Display ..... \$15.00
- RFM Wireless (RF) Mouse Only Module ..... \$750.00  
**RFM will be available in September, 2010. Allows wireless computer mouse operation through the driver control. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ joysticks.**
- IR Infrared Module with Wireless RF Mouse Control ..... \$995.00  
**Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.**

## ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black ..... \$156.00
- ASL611 Adjustable Switch Mounting Bracket ..... \$252.00

## MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks (Coming Soon) ..... \$45.00
- VAN02 Oxygen Holder (Van Seat Only) ..... \$150.00 E2208

## BATTERY TRAY TYPE

- 24TRY 24 Group Style Tray ..... No Charge

## BATTERY

- 24BATTERY 24 Gel Battery (On Chair) ..... \$410.00 E2363

## BATTERY CHARGER OPTION

- Clear
- 110CHARGER 110 Volt Battery Charger ..... No Charge
- LESTR Lester 8 AMP Battery Charger ..... \$320.00 E2367
- L8900 Omit Battery Charger ..... Less \$80.00

## CENTER MOUNT

Center mount adjustable angle (14"-18") available with tall STF only. 13.25" - 17.25" extension range available through special #SM-071217095848

- CTMT-F Center Mount Adj. Knee Angle (11" - 15") ..... \$300.00  
**Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 11" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.**
  - SM Small Footplates (5 1/4" x 7 1/2") ..... No Charge
  - MED Medium Footplates (5 1/2" x 9 1/2") ..... No Charge
  - LRG Large Footplates (5.90" x 11 1/2") ..... No Charge
  - CALF Calf Pads For Center Mount ..... \$200.00
- CTMT-FJ Jr. Center Mount Adj. Knee Angle (7 1/2" - 11") ..... \$300.00  
**Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 7 1/2" - 11". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.**
  - SM Small Footplates (5 1/4" x 7 1/2") ..... No Charge
  - MED Medium Footplates (5 1/2" x 9 1/2") ..... No Charge
  - LRG Large Footplates (5.90" x 11 1/2") ..... No Charge
  - CALF Calf Pads For Center Mount ..... \$200.00

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® TDX® SR Wheelchair with Captain's Seat

## CENTER PIVOT STYLE (PAIR)

Clear

All front riggings include a standard footplate.

- 70NHD 70° Swingaway Footrests 13 1/2" - 17 1/2" ..... Standard  
With seat widths of 21" to 24" and depth of 16", must pick footplate 1350.
- 70STAPER 70° Tapered Footrests 13" - 17" ..... No Charge  
Standard with impact guards.
- AT5544 Manual Elevating Legrests ..... \$255.00 E0990  
Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

## PIN STYLE-HEMI (PAIR)

Clear

All front riggings include a standard footplate.

- PHW93 70° Swingaway Footrests 13 1/2" - 19" ..... No Charge
- PHAL4A Manual Art./ELR 17 1/4" - 21 1/2" ..... \$225.00 K0053

## FOOTBOARD

- 90FB Footboard ..... \$150.00

## FOOTPLATE OPTIONS (PAIR)

Clear

- 1651 Composite Flip-Up Footplates ..... No Charge
- 1350 Extra Large Aluminum Footplates ..... No Charge  
Not available with AT5544 and 17" wide by 16" deep and low seat-to-floor.
- AT5543 Adjustable Angle Flip-Up Footplates ..... \$165.00 K0040  
Not available with AT5544 and 16" wide by 16"-18" deep or 17" wide by 16" deep and low seat-to-floor.
- CTMTFB Footboard for Center Mount ..... No Charge  
One piece footboard for Center Mount (CTMT-F) Rigging.

## FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap ..... \$27.00 E0038
- IMPACTG Impact Guards ..... \$51.00
- HLOOP Heel Loops (pair) ..... \$45.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 ..... \$83.00
- ALPT Longer Pivot And Slide Tube ..... \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... \$84.00

TDXSR Specifications	
Seat Width Range	Adult (Adjustable in 1" increments) 16", 18", 20", 22"
Seat Depth Range	16" - 22"
Back Height Range	12" - 24"
Back Angle Range	85° to 105° in 5" increments
Seat-to-Floor (approx)	21" Standard
Overall Width of Base	24" (without Joystick)
Overall Height	(with 20" back height and 0° seat angle)
Semi-Recline w/o headrest	38"
Semi-Recline w/ headrest	45"
Overall Length	
without front rigging	35 1/4"
with 90° front rigging	43"
Drive Axle	Non-Adjustable
Drive Wheels/Tires (foam filled)	
Standard	14" x 3"

### Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

### Weight Limit 300 or 400 Pounds

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**Note:** All weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

Comments

**NOTE: All specifications and dimensions are approximate.**  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).



**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)

USA

One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

©2013 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.