

I N V A C A R E

POWER CHAIRS

A **Quick Reference Guide** to Invacare Power Chairs



REAR Wheel Drive
AVIVA® Storm RX™



MID Wheel Drive
ROVI X3 and A3



MID Wheel Drive
TDX® SP2



FRONT Wheel Drive
AVIVA® FX





Pelvic Positioning Strap not shown on chairs for clarity purposes.

Design. Technology. Performance.

At Invacare Corporation we make Power Wheelchairs with Seating that help improve the lives of our customers. By listening to our customers and their care teams, we focus on meaningful design and purposeful technology that enhances performance.



POWER CHAIRS

Chair Base/Drive Type	AVIVA® Storm RX™ REAR Wheel Drive	ROVI X3 MID Wheel Drive	ROVI A3 MID Wheel Drive	TDX® SP2 MID Wheel Drive	AVIVA® FX FRONT Wheel Drive
					
Base Weight Capacity	300 lbs. - Standard 250 lbs. - Narrow 175 lbs. - Ultra Low Mini	300 lbs.	300 lbs.	300lbs. 400lbs. HD w/ Ultra Low Maxx, 450lbs w/captain seat	300 lbs.
HCPCS Codes	K0848, K0856, K0861	K0848, K0849, K0856, K0861	K0848, K0849, K0856, K0861	K0848, K0849, K0850, K0851, K0856, K0858, K0861, K0862	K0848, K0849, K0856, K0861
Speed	Standard - 5.8 mph Performance - 7.5 mph (Std width only)	Standard - 5.0 & 5.8 mph Performance - 6.2 mph	Standard - 5.8 & mph Performance - 6.2 mph	Standard - 5.8 mph Performance - 7.5 mph	Standard - 5.8 mph Performance - 7.5 mph
Drive Wheels	14"	14"	14"	14"	14"
Caster Size	8" x 2.5" Black Casters w/Solid Inlay	6" Solid (Front/Rear)	6" Solid (Front/Rear)	6" Solid (Front/Rear)	8" x 2.5" Black Casters w/Solid Inlay
Ground Clearance	>3.2" (under battery box)	3"	3"	2.5"	3.15"
Minimum Base Turning Radius	24.8" - Narrow 25" - Standard	19.5"	20.5"	20"	26.8"
Driving Range	Up to 19.7 miles - Standard Up to 13.8 miles - Narrow	Up to 17.7 miles	Up to 17.7 miles	Up to 20 miles	Up to 20 miles
Base Width	22.4" - Narrow 24.4" - Standard	23.25"	25.5"	24" & 25.5"	24.25"
Base Length	33.5"	36.5"	36.6"	35.6"	42.1"
Suspension	SureGrip™ Suspension System	Active Ride Control (ARC)	Active Ride Control (ARC)	Surestep® and Stability Lock Suspension	4Sure™ Suspension
Tracking Technology	G-Trac™ - Digital Gyroscopic Technology (Optional)	OnTraxx - Enhanced Driving Module	OnTraxx - Enhanced Driving Module	G-Trac - Digital Gyroscopic Technology (Optional)	G-Trac - Digital Gyroscopic Technology
Motor Type	DuraWatt™ 4-Pole Sealed Housing	4-Pole EAD Motors	4-Pole EAD Motors	DuraWatt 4-Pole Sealed Housing	DuraWatt 4-Pole Sealed Housing
Battery Type	Group 24 - Standard M50 - Narrow	Group 34 & Group 50	Group 34 & Group 50	22NF & Group 24	Group 24
Battery Isolation Switch	YES	NO	NO	NO	YES
Incline Climbing	9 degrees	8 degrees	8 degrees	9 degrees	9 degrees
EZ Lock Compatible	YES	YES	YES	YES	YES
Lights	LED integrated into the front caster and a rear LED strip built into the rear shroud. (standard)	Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. (optional)	Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. (optional)	LED Light Package Available for front and rear lights. (optional)	LED Light package available for front and rear lights. (optional)

Chair Base/Drive Type	AVIVA® Storm RX™ REAR Wheel Drive	ROVI X3 MID Wheel Drive	ROVI A3 MID Wheel Drive	TDX® SP2 MID Wheel Drive	AVIVA® FX FRONT Wheel Drive
-----------------------	-----------------------------------	-------------------------	-------------------------	--------------------------	-----------------------------

SEATING					
Captain Seat	NO	YES	YES	YES	YES
Rehab Seating	YES	YES	YES	YES	YES
Ultra Low Mini Maxx Power Positioning System	YES	YES	NO	YES	YES
Ultra Low Maxx Power Positioning System	YES	YES	YES	YES	YES
Ultra Low Bariatric Power Positioning System	NO	NO	NO	YES/HD	NO
MPS Mini Maxx - Power Standing System	NO	YES	YES	NO	NO
MPS Maxx - Modular Power Standing System	NO	YES	YES	NO	YES
407 Maxx - Single Pivot Tilt Power Positioning System	NO	YES	YES	NO	NO

ELECTRONICS					
Technology	LiNX®	R-Net	R-Net	LiNX®	LiNX®
Driver Control Options	<ul style="list-style-type: none"> ▶ LED Joystick ▶ 3.5" Color Touch Screen ▶ Color Touchscreen Display 	<ul style="list-style-type: none"> ▶ LED Joystick ▶ Color Remote Joystick ▶ Advanced Joystick ▶ Omni Display 	<ul style="list-style-type: none"> ▶ LED Joystick ▶ Color Remote Joystick ▶ Advanced Joystick ▶ Omni Display 	<ul style="list-style-type: none"> ▶ LED Joystick ▶ 3.5" Color TouchScreen ▶ Color Touchscreen Display 	<ul style="list-style-type: none"> ▶ LED Joystick ▶ 3.5" Color Touch Screen ▶ Color Touchscreen Display
Remote Diagnostics	MyLiNX	N/A	N/A	MyLiNX	MyLiNX
Bluetooth	YES	YES	YES	YES	YES

OTHER FEATURES					
Limited Warranty	<ul style="list-style-type: none"> ▶ 1 Year Electronics/ Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries 	<ul style="list-style-type: none"> ▶ 1 Year Electronics/ Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries 	<ul style="list-style-type: none"> ▶ 1 Year Electronics/ Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries 	<ul style="list-style-type: none"> ▶ 1 Year Electronics/ Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries 	<ul style="list-style-type: none"> ▶ 1 Year Electronics/ Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries

COLORS														
Obsidian Black	Ruby Red	Denim Navy	Black Carbon (Matte)	Mellow Yellow	Motion Green	Black Carbon (Matte)	Mellow Yellow	Motion Green	Black Ice (Gloss)	Space Station Silver	Invacare Blue	Obsidian Black	Ruby Red	Denim Navy
Morning Mist	Orange Marmalade	Pink Flamingo	Orange Crush	Indigo Blue	Party Pink (Metallic)	Orange Crush	Indigo Blue	Party Pink (Metallic)	Rockstar Red	Island Blue	Tangy Orange	Morning Mist	Orange Marmalade	Pink Flamingo
Lime Twist	Purple Passion	Arctic Blue	Nordic Blue	Purple Reign	Blizzard White (Metallic)	Nordic Blue	Purple Reign	Blizzard White (Metallic)	Grape Jelly Bean	Lights Out Black (Matte)	Monster Green	Lime Twist	Purple Passion	Arctic Blue
Blue Peacock	Banana Split		Candy Apple Red (Metallic)	Graphite Gray	Black Panther (Gloss)	Candy Apple Red (Metallic)	Graphite Gray	Black Panther (Gloss)	White Out			Blue Peacock	Banana Split	

POWER CHAIR ADD-ONS



Invacare® Matrix® Seating & Positioning

From our highly regarded Libra™ cushion and award winning E2™ backs, our Mini pediatric line to Matrix bariatric products, we offer simple, safe and effective seating solutions to help promote the quality of life of our customers.



PinDot Custom Molded Seating

Custom molded seating solutions for end users with complex seating needs. State-of-the-art technology, and computer assisted manufacturing helps us create complex seating for each and every customer.



Adaptive Switch Labs

ASL designs products specifically for those who cannot access their environment through conventional means. Our goal has always been independence and accessibility for those with whom we work.



[REHAB.INVACARE.COM](https://rehab.invacare.com)

©2022 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. LINX is a registered trademark of Dynamic Controls. Specifications are subject to change without notification. Form no. 21-154 rev 4.22 220417 Please refer to the User's Manual for additional safety information and instructions.

