



ALBER POWER ASSIST

CHOOSE YOUR POWER



MOBILITY LEVEL
ACTIVE

SMOOV[®] one

Innovative electric drive that easily docks on and off your wheelchair for when you need power.

Aluminum locking claw helps the unit stay secure to the chair – no “pop off”

Rubber tire works on all types of terrain and surfaces including gravel, wet roads or cobblestones

356 degree swiveling fork Provides power into and while in turns for enhanced performance

Speed up, slow down and stop easily with the ergonomic control unit

App allows you to customize your ride, check battery life, record trips and more...

Includes the drive unit, control unit, chair connection, charger, and app.

SPECIFICATIONS

- Max User Weight: 310 lbs.
- Max User Weight Plus Chair & Accessories: 374 lbs.
- Max Speed: 6.0 mph*
- Max Incline: Up to 16%
- Range: Up to 12.3 Miles*
- Weight of Drive Unit: 16 lbs.



MOBILITY LEVEL
MEDIUM ACTIVE

twion[®]

Wheels with in-hub motor and battery, that propel the chair in proportion to the user’s push.

Easy to transport - detach the wheels with our quick release feature

Features high-tensile aluminum alloys, high-tech plastics and the latest digital motor technologies

Lithium ion battery – charges up to 80% in 4 hours

Provides proportional power – the harder a user pushes, the more power they receive

Includes two wheels with built in motor & lithium ion battery, and charger.

SPECIFICATIONS

- Max User Weight: 250 lbs.
- Max User Weight Plus Chair & Accessories: 320 lbs.
- Max Speed: 6.0 mph*
- Max Incline: Up to 16%
- Range: Up to 12 Miles*
- Weight of Each Wheel: 13 lbs.
- Wheel Size: 24”



MOBILITY LEVEL
MEDIUM ACTIVE

e-motion[®]

Wheels with in-hub motor and battery, that provide the user with additional power.

Easy to transport - detach the wheels with our quick release feature

The newly designed motor and battery create more power and efficiency

Includes cruise mode which keeps a constant speed with a single push – just use your hands to steer the chair. Works on both inclines and declines (with optional app upgrade)

Patented roll-back delay holds position for 5 seconds after release of the hand rim while driving uphill

Assisted braking while going downhill

App allows you to customize your ride, check battery life, record trips and more...

Includes two wheels with built in motor & lithium ion battery, charger, app, and optional control unit (if you are not using the app)

SPECIFICATIONS

- Max User Weight: 330 lbs.
- Max User Weight Plus Chair & Accessories: 418 lbs.
- Max Speed: 5.3 mph*
- Range: Up to 15.5 Miles*
- Weight of Each Wheel: 17 lbs.
- Wheel Size: 22”, 24”, 25”



MOBILITY LEVEL
LOW ACTIVE

e-fix[®]

Add-on drive with joystick that helps turn your manual wheelchair into a portable power wheelchair

Program your e-fix to meet your needs by using the compact control unit with swivel arm and color display

Portable power – each of the lightweight components can be removed quickly and easily without tools

Powerful drive wheels work great indoors and outdoors

Attendant control option allows a care giver to control the power and braking

Cover longer distances quickly and easily

Includes wheels with motor, battery pack, joystick with control display, charger, and optional attendant control unit

SPECIFICATIONS

- Max User Weight: 350 lbs.
- Max User Weight Plus Chair & Accessories: 460 lbs.
- Max Speed: 3.75 mph*
- Range: Up to 10 miles (standard) – 12.4 miles (with add on high capacity battery)*
- Weight of each drive wheel: 17 lbs.
- Weight of all components: 42 lbs.
- Wheel Size: 24”

Please refer to the Alber and wheelchair user manuals for additional safety information and instructions.

*Varies based on weight of the user, temperature and terrain

CONTACT YOUR INVACARE SALES REPRESENTATIVE TO LEARN MORE OR SCHEDULE A DEMO.
888.426.8581 | info@alber-usa.com | alber-usa.com | Alber USA LLC | 1005 International Drive | Oakdale, PA 15071

